

The background of the slide features a close-up, shallow depth-of-field photograph of a multi-well plate. A pipette tip is positioned above one of the wells, dispensing a drop of bright pink liquid. The other wells in the plate also contain a similar pink liquid. The lighting is soft and focused on the pipette and the liquid, creating a clean, scientific aesthetic.

PacBio

# Technical overview: Whole genome and metagenome library preparation using SMRTbell prep kit 3.0

Sequel II and IIe Systems ICS v11.0 / SMRT Link v11.0

PN 102-390-900 Version 01 (April 2022)

# Preparing whole genome and metagenome libraries using SMRTbell prep kit 3.0

## Technical overview

1. WGS library sample preparation workflow overview
2. WGS library sample preparation workflow details
3. WGS library sequencing preparation workflow overview
4. WGS data analysis recommendations for variant detection and *de novo* assembly applications
5. WGS library example sequencing performance data
6. Technical documentation & applications support resources
7. APPENDIX – Genomic DNA isolation & QC recommendations for PacBio WGS library sample preparation

# Whole genome sequencing: How to get started

Application-specific  
Best practices guide

Application-specific  
Procedure & checklist

Application-specific  
technical overviews

Library construction,  
sequencing & analysis

**WHOLE GENOME SEQUENCING FOR *DE NOVO* ASSEMBLY – BEST PRACTICES**

PacBio® HiFi reads provide both long read lengths (up to 25 kb) and high accuracy (>99.9%) to quickly and affordably generate contiguous, complete, and correct *de novo* genome assemblies of even the most complex genomes.

**Contiguity**  
High contig N50

**Completeness**  
No missing bases or fragmented genes

**Correctness**  
High base accuracy and phased alleles

**Compute**  
Small file sizes and fast analysis time

The Sequel® IIe system provides cost-effective and scalable HiFi sequencing of any genome

**Large or complex genomes**

Fast and efficient assembly of even the largest genomes with haplotype resolution of complex polyploids

The 27 Gb hexaploid genome of the redwood tree was sequenced and assembled in under two weeks

**Human genomes at scale**

Flexible and scalable workflows for sequencing 100s–1000s of human genomes per year from a variety of sample types

Assemble a human genome in one day

"If your genome isn't HiFi, it's no longer reference grade."  
Kevin McKernan, Medicinal Genomics

**Preparing whole genome and metagenome libraries using SMRTbell® prep kit 3.0**

Procedure & checklist

**Before you begin**

This procedure describes the workflow for constructing whole-genome sequencing (WGS) libraries from genomic and metagenomic DNA using the SMRTbell prep kit 3.0 for sequencing on PacBio systems.

Overview	
Samples per SMRTbell prep kit 3.0	1–24
Workflow time	4.5 hours for up to 8 samples; 6 hours for 24 samples Time difference is from DNA shearing, which is done in sets of 8 samples. Excludes measuring DNA size on Femto Pulse system.

DNA Input		Microbes	Metagenomes
Quantity	300 ng–5 µg per library		
DNA size distribution (Femto Pulse system)	50% ≥ 30 kb & 90% ≥ 10 kb	90% ≥ 7 kb	90% ≥ 7 kb
DNA Shearing (MagneSpirator 3 system)	Speed 31	Speed 40	Speed 40
Target fragment lengths	15–18 kb	7–12 kb	7–12 kb
Size selection required	AMPure® PB beads	none	none

© 2022 PacBio. All rights reserved. Research use only. Not for use in diagnostic procedures.  
PN 102-166-600. EA11 18FEE2022

**Technical Overview: Preparing whole genome and metagenome libraries using SMRTbell prep kit 3.0**

Sequel II and IIe Systems ICS v11.0 / SMRT Link v11.0

PN 102-390-900 Version 01 (April 2022)

**Example library QC and sequencing performance for human WGS libraries prepared with SMRTbell prep kit 3.0**

SMRTbell library QC and primary sequencing metrics

**Size-Selected Library QC**

1x 1 kb Final Library (Mean Size Range: 1.0 kb, Min: 0.8 kb, Max: 1.2 kb)

Identify size for 15–18 kb target (blue line) and target (red line)

**Raw Data Report**

80 kb  
10 kb

80% Read Length Accuracy

**CCS Analysis Report**

14.8 kb Mean HiFi Read Length

99.99% Accuracy

Input gDNA for MagneSpirator 3 shearing	3000 ng	Raw Base Yield	637.65 Gb	HiFi Reads	2.7 M
Post-shearing recovery (%)	2620 ng (87%)	Mean Polymerase Read Length	102.8 kb	HiFi Base Yield	39.2 Gb
Final yield of AMPure PB bead size-selected library (%)	1070 ng (36%)	P1	29.4%	Mean HiFi Read Length	14,460 bp
		P2	75.0%	Median HiFi Read Quality	Q34
		P3	1.4%	HiFi Read Mean # of Passes	12

© 2022 PacBio. All rights reserved. Research use only. Not for use in diagnostic procedures.  
PN 102-166-600. EA11 18FEE2022

Application Briefs: [WGS for \*de novo\* assembly – Best practices](#) / [Variant detection using WGS with HiFi Reads – Best practices](#)


Summary overview of application-specific sample preparation and data analysis workflow recommendations

Procedure & Checklist – Preparing whole genome and metagenome libraries using SMRTbell prep kit 3.0 (102-166-600)

Technical documentation containing sample library construction and sequencing preparation protocol details

Technical overview: Whole genome and metagenome library preparation using SMRTbell prep kit 3.0 (102-390-900)

Technical overview presentations describe sample preparation details for constructing HiFi libraries for specific applications. Example sequencing performance data for a given application are also summarized.



**Genomic DNA QC & shearing**

≥1 µg per SMRT Cell 8M  
Shear to 15 – 18 kb for large genomes  
Shear to 10 – 12 kb for small genomes

**Library construction (SMRTbell Prep Kit 3.0)**

Multiplex WGS samples using SMRTbell barcoded adapter plate 3.0

**Sequencing (Sequel II and IIe systems)**

ABC\* with Sequel II binding kit 3.2  
15 hr or 30 hr movie collection time

**Data analysis (SMRT Link or third-party)**

Genome assembly  
Variant detection

# When is it appropriate to consider standard input, low DNA input or ultra-low DNA input workflows for whole genome sequencing applications?

Sample & project considerations	Standard HiFi WGS (Large/small genomes & shotgun metagenomics)	Low DNA input WGS (2-plex)	Low DNA input WGS (Single sample)	Ultra-low DNA input WGS
<b>Recommended Procedure &amp; checklist</b>	Preparing whole genome and metagenome libraires using SMRTbell prep kit 3.0 <a href="tel:102-166-600">(102-166-600)</a>	Preparing whole genome and metagenome libraires using SMRTbell prep kit 3.0 <a href="tel:102-166-600">(102-166-600)</a>	Preparing whole genome and metagenome libraires using SMRTbell prep kit 3.0 <a href="tel:102-166-600">(102-166-600)</a>	Preparing HiFi SMRTbell Libraries from Ultra-Low DNA Input <a href="tel:101-987-800">(101-987-800)</a>
<b>Minimum DNA input</b>	≥3 µg for a 3-Gb genome	300 ng for each genome	400 ng	5 ng
<b>Amplification-based?</b>	No	No	No	Yes
<b>Genome size limit</b>	N/A	600 Mb for each genome	1 Gb	500 Mb
<b>Supported applications</b>	<i>De novo</i> Assembly Human Variant Detection	<i>De novo</i> Assembly	<i>De novo</i> Assembly	<i>De novo</i> Assembly Human Variant Detection

**Ultra-low DNA input:**  
**SUPPORTED APPLICATIONS**



**Ultra-low DNA input:**  
**UNSUPPORTED APPLICATIONS**






De novo assembly of insect/arthropod genomes (Up to 500 Mb)




Variant detection (SNPs, Indels, SVs) in human genomes (3 Gb)

De novo assembly for microbes, plants, vertebrates, or other non-DNA limited sample types

Metagenomics sequencing



# WGS library sample preparation workflow overview

# WGS sample preparation procedure description

Procedure & Checklist – Preparing whole genome and metagenome libraries using SMRTbell prep kit 3.0 ([102-166-600](https://www.pacb.com/support/102-166-600)) describes a method for constructing SMRTbell libraries that are suitable for generating HiFi reads on the Sequel II and IIe Systems for **WGS and metagenomic shotgun sequencing applications**.

## Procedure Highlights

- Uses **SMRTbell Prep Kit 3.0** (102-182-70) and supports high-throughput processing using **500 ng – 5 µg** of input genomic DNA amounts
  - We recommend starting with **≥1 µg of input DNA per SMRT Cell 8M** (or ~3 µg for up to a 3 Gb WGS sample to enable running 3 SMRT Cells 8M)
- Multiplexing of samples can be performed using **SMRTbell barcoded adapter plate 3.0** (102-009-200)
- Recommend shearing high-quality gDNA using a **Megaruptor 3 System** (Diagenode)
  - **15 kb – 18 kb** target insert size for large (plant / animal / human) genomes
  - **7 kb – 12 kb** target insert size for small (microbial) genomes
  - **7 kb – 12 kb** target insert size for shotgun metagenomic samples
- **4.5-hour workflow time** to process up to 8 samples from shearing to size selection (6 hours for 24 samples)
  - Time difference is from DNA shearing, which can be performed in sets of 8 samples.
  - Excludes time needed for DNA sizing QC analysis using a Femto Pulse system.
- WGS SMRTbell libraries can be **size-selected using AMPure PB Beads** without the need for third-party equipment

Preparing whole genome and metagenome libraries using SMRTbell® prep kit 3.0

PacBio

Procedure & checklist

## Before you begin

This procedure describes the workflow for constructing whole-genome sequencing (WGS) libraries from genomic and metagenomic DNA using the SMRTbell prep kit 3.0 for sequencing on PacBio systems.

Overview			
Samples per SMRTbell prep kit 3.0	1–24		
Workflow time	4.5 hours for up to 8 samples; 6 hours for 24 samples Time difference is from DNA shearing, which is done in sets of 8 samples. Excludes measuring DNA size on Femto Pulse system.		
DNA Input			
Quantity	300 ng–5 µg per library		
	Human, plant, and animal	Microbes	Metagenomes
DNA size distribution (Femto Pulse system)	50% ≥ 30 kb & 90% ≥ 10 kb	90% ≥ 7 kb	90% ≥ 7 kb
DNA Shearing (Megaruptor 3 system)	Speed 31	Speed 40	Speed 40
Target fragment lengths	15–18 kb	7–12 kb	7–12 kb
Size selection required	AMPure® PB beads	none	none

© 2022 PacBio. All rights reserved. Research use only. Not for use in diagnostic procedures. PN 102-166-600 EA V1 18FEB2022

PacBio

PacBio [Documentation](https://www.pacb.com/support/102-166-600) (102-166-600)


APPLICATIONS

WHOLE GENOME SEQUENCING

*De Novo assembly & variant detection*

*Microbial assembly*

*Shotgun metagenomics*



# WGS sample preparation & sequencing workflow overview

Workflow summary for constructing SMRTbell libraries suitable for sequencing on the Sequel II and Ie systems for WGS and metagenomic shotgun sequencing applications




## Genomic DNA extraction, QC & shearing

- Perform DNA sample extraction using [Circulomics Nanobind kits](#)
- $\geq 1$   $\mu\text{g}$  of total input DNA per SMRT Cell 8M (for a single sample or across multiple samples when pooling)
- Large genome samples: Femto Pulse Genome Quality Number (GQN) at 30 kb  $\geq 5.0$
- Small genome samples: GQN at 7 kb  $\geq 9.0$
- Shear to target insert size with Megaruptor 3 system



## Sequencing

- Follow **SMRT Link Sample Setup** instructions for primer annealing, polymerase binding, complex cleanup and sample loading



HiFi Read

PacBio HiFi reads achieve 99.9% accuracy




## SMRTbell library construction

- **Procedure & checklist – Preparing whole genome and metagenome libraries using SMRTbell prep kit 3.0** ([102-166-600](#))
- Multiplex WGS samples using SMRTbell barcoded adapter plate 3.0 ([102-009-200](#)).
- Perform AMPure PB bead size selection for large genome WGS libraries or perform a standard cleanup using SMRTbell cleanup beads for small genome WGS libraries



## Data analysis

- For *de novo* assembly, can use [SMRT Link Genome Assembly](#) or other third-party software
- For variant detection, can use [DeepVariant](#) for small variants  $< 20$  bp and SMRT Link PBSV for larger variants  $> 20$  bp





# WGS library sample preparation workflow details



# Procedure & checklist – Preparing whole genome and metagenome libraries using SMRTbell prep kit 3.0

Procedure & checklist [102-166-600](#) describes a method for constructing SMRTbell libraries using SMRTbell prep kit 3.0 that are suitable for generating high-accuracy long reads on the Sequel II and IIS systems for **whole genome sequencing** (*de novo* assembly, variant detection, microbial multiplexing) and **shotgun metagenomic sequencing**

## Protocol Contents

1. Recommendations for **gDNA quantification and sizing QC**.
2. Recommendations for shearing gDNA to the desired target mode size using the **Megaruptor 3** system (Diagenode).
3. Enzymatic steps for preparation of a WGS SMRTbell library using **SMRTbell prep Kit 3.0** (102-182-700). (Instructions for preparing **multiplexed** samples using **SMRTbell barcoded adapter plate 3.0** (102-009-200) are also provided.)
4. Instructions for size-selection of WGS SMRTbell libraries using **AMPure PB bead size selection**. (Size selection is not required for microbial WGS and metagenomic shotgun libraries where retention of shorter fragments is desired.)
5. Guidance for **pooling** barcoded WGS SMRTbell libraries for multiplexed sequencing on a single SMRT Cell.

## Preparing whole genome and metagenome libraries using SMRTbell® prep kit 3.0

PacBio

Procedure & checklist

### Before you begin

This procedure describes the workflow for constructing whole-genome sequencing (WGS) libraries from genomic and metagenomic DNA using the SMRTbell prep kit 3.0 for sequencing on PacBio systems.

Overview			
Samples per SMRTbell prep kit 3.0	1–24		
Workflow time	4.5 hours for up to 8 samples; 6 hours for 24 samples Time difference is from DNA shearing, which is done in sets of 8 samples. Excludes measuring DNA size on Femto Pulse system.		
DNA input			
Quantity	300 ng–5 µg per library		
	Human, plant, and animal	Microbes	Metagenomes
DNA size distribution (Femto Pulse system)	50% ≥ 30 kb & 90% ≥ 10 kb	90% ≥ 7 kb	90% ≥ 7 kb
DNA Shearing (Megaruptor 3 system)	Speed 31	Speed 40	Speed 40
Target fragment lengths	15–18 kb	7–12 kb	7–12 kb
Size selection required	AMPure® PB beads	none	none

© 2022 PacBio. All rights reserved. Research use only. Not for use in diagnostic procedures.  
PN 102-166-600 EA V1 18FEB2022

PacBio

PacBio [Documentation](#) (102-166-600)

# SPK 3.0 WGS library prep and sequencing workflow timing overview


Go from DNA shearing to sequencing preparation in 1 day with SPK 3.0

Workflow step		Hands-on (min)	Walk-away (min)
DNA shearing*	DNA shearing (1 Megaruptor 3 cycle)	5	45
	1.0X SMRTbell bead cleanup	5	20
SMRTbell library construction*	DNA repair / a-tailing	5	35
	Adapter ligation (barcoded or non-barcoded adapter)	5	30
	1.0X SMRTbell bead cleanup	5	20
	Nuclease treatment	5	15
	AMPure PB bead size selection	10	30
<b>Total</b>		<b>40 min</b>	<b>3.3 hrs</b>

Workflow step		Hands-on (min)	Walk-away (min)
Sequencing preparation (ABC)	Primer annealing (Sequel II primer 3.2)	5	15
	Polymerase binding (Sequel II binding kit 3.2)	5	15
	Complex cleanup (1.2X SMRTbell cleanup beads)	5	20
<b>Total</b>		<b>15 min</b>	<b>0.83 hrs</b>

## SPK 3.0 WGS Workflow

(102-166-600)



# Input genomic DNA QC recommendations for WGS library construction using SMRTbell Prep Kit 3.0

- WGS library construction using SMRTbell prep kit 3.0 requires **high-quality, high-molecular weight genomic DNA\***.
- Prior to library preparation, evaluate the **quantity** and **size distribution** of the input gDNA to determine whether it is suitable for the protocol.
- For each input gDNA sample:
  - Measure concentration and total mass of DNA with a **Qubit High Sensitivity dsDNA Assay** system (Thermo Fisher Scientific)
  - Measure DNA size distribution with a **Femto Pulse** system (Agilent)
  - Proceed with SMRTbell library construction if the **gDNA sample quality** is acceptable (see Table below)

Sample type	Input DNA metric	Requirement	Notes
All	Per Library	300 ng – 5 µg	<ul style="list-style-type: none"> <li>• Starting with low DNA input amounts approaching <b>~300 ng</b> may in some cases produce insufficient amounts of SMRTbell library to load at concentrations that optimize sequencing data yield.</li> <li>• For multiplexing applications, generally aim to <b>use ≥300 ng of DNA input per sample</b>, with a total mass ≥1 µg across all samples</li> </ul>
All	Per SMRT Cell 8M	≥1 µg	<ul style="list-style-type: none"> <li>• Start with <b>≥1 µg of total input DNA per SMRT Cell 8M</b> (for a single sample or across multiple samples when pooling) to enable generation of sufficient library to load at concentrations that optimize sequencing data yield.</li> </ul>
Large genome (Animal/plant/human)	Longer than 30 kb	≥50%	<ul style="list-style-type: none"> <li>• Required to achieve target fragment lengths after DNA shearing.</li> <li>• <b>For large genome samples, the Femto Pulse Genome Quality Number (GQN) at 30 kb should be ≥5.0.</b> (Not applicable to microbial and metagenomic samples)</li> </ul>
Large genome (Animal/plant/human)	Longer than 10 kb	≥90%	<ul style="list-style-type: none"> <li>• Required for effective AMPure PB bead size selection.</li> <li>• <b>For large genome samples, the GQN at 10 kb should be ≥9.0.</b></li> </ul>
Small genome (microbial/metagenomic)	Longer than 7 kb	≥90%	<ul style="list-style-type: none"> <li>• <b>For microbial and metagenomic samples, the input DNA should be at least as large as the recommended insert lengths of 7–12 kb with a GQN at 7 kb ≥9.0.</b></li> <li>• Any degradation present should be due to shearing from the extraction process (e.g., bead beating) and not from poor sample handling or storage, or biochemical processes</li> </ul>

# DNA extraction tech note: Sample preparation for PacBio HiFi sequencing from human whole blood (102-326-500)


Provides best practices for handling human whole blood samples\* to generate optimal sequencing performance

Technical note [102-326-500](#) Discusses the effect of anticoagulant, sample storage time, storage conditions, and white blood cell count on the sequencing performance of DNA extracted using the Nanobind CBB Big DNA Kit ([NB-900-001-01](#))



Nanobind CBB Big DNA Kit (NB-900-001-01) for isolating HMW DNA from cells, bacteria, & blood.

For optimal HiFi yield and read length performance, store human whole blood samples for fewer than 2 days at 4°C.




Stage	Variable	Best practice for PacBio HiFi sequencing
Before DNA extraction	Sample type	Human whole blood
	Anticoagulant	Potassium EDTA (K <sub>2</sub> EDTA)
	Sample storage temperature	4 ± 3°C
	Sample storage time	≤ 2 days from collection to extraction
DNA extraction	Volume of whole blood	200 µL
	White blood cell (WBC) count	≥ 4 × 10 <sup>6</sup> cells/mL for ≥ 3 µg of DNA
	DNA extraction method	Nanobind CBB Big DNA kit
After DNA extraction	DNA storage	Rest 1 day at ambient temperature, then store at 4 ± 3°C
	DNA size distribution	<ul style="list-style-type: none"> <li>90% of DNA ≥ 10 kb (genomic quality number at 10 kb ≥ 9.0)</li> <li>50% of DNA ≥ 30 kb (genomic quality number at 30 kb ≥ 5.0)</li> </ul>
	UV absorbance	<ul style="list-style-type: none"> <li>A<sub>260/280</sub> nm ≥ 1.7</li> <li>A<sub>260/230</sub> nm ≥ 1.5</li> </ul>

# The Femto Pulse System is recommended for DNA sizing QC of genomic DNA for WGS applications

- Femto Pulse System (Agilent) is highly recommended for DNA sizing QC of input genomic DNA and SMRTbell libraries
  - Enables sizing of gDNA samples ranging from 1,300 bp to 165 kb
  - Requires <1 ng of sample DNA
  - Can analyze up to 12 samples in <1.5 hours
- The Femto Pulse system can be used in place of traditional pulse-field gel electrophoresis (PFGE) to quickly assess the initial integrity of genomic DNA, evaluate shears, determine appropriate size-selection thresholds, and conduct final QC before preparing libraries for SMRT Sequencing

Femto Pulse outputs quality metrics such as the **Genomic Quality Number (GQN)\*** to quickly score the integrity of HMW gDNA

Ave. Smear Size (bp)	GQN Set at 30 kb
12,147	0
23,339	1.5
45,304	6.4
57,789	7.1
73,267	7.8
94,045	7.8
109,968	8.2
164,292	8.8



Femto Pulse System



Femto Pulse system offers a simplified QC workflow to generate SMRTbell libraries for WGS sequencing in **reduced time**, **and conserves sample** by using femtogram ranges of input DNA

# The Megaruptor 3 system is recommended for shearing genomic DNA For WGS applications

- Megaruptor 3 system (Diagenode) is highly recommended for DNA shearing\*
  - Up to 8 samples can be sheared in parallel in ~45 minutes for high-throughput applications
  - Achieving the same size distribution across multiple samples provides more consistent sequencing performance
- Recommended library insert size distributions and Megaruptor 3 shear speed settings to use for different WGS applications are summarized on Page 7 in the procedure
  - Bring input gDNA to a final volume of **100 – 130  $\mu\text{L}$  with Low TE buffer** [10 mM Tris-HCl (pH 8.0) + 0.1 mM EDTA] to target a DNA concentration of **3 – 39 ng/ $\mu\text{L}$  (ideal: 30 ng/ $\mu\text{L}$ )**
  - Perform shearing (**1 cycle**) using the conditions described in the table below

Megaruptor 3 System



Application	Recommended Library Insert Size (Mode)	Recommended Megaruptor 3 Shear Speed Setting
Animal / plant / human WGS	15 kb - 18 kb	31
Microbial WGS or shotgun metagenomics	7 kb - 12 kb	40


Because the response of individual gDNA samples can differ, **optimization of shearing conditions** may be needed to achieve the desired fragment distribution



\* **Note:** The g-TUBE (Covaris) device generates a broader DNA fragment size-distribution compared to the Megaruptor 3 system. As a result, HiFi read quality and overall HiFi data yield may be reduced due to the residual presence of very large DNA fragments generated by g-TUBEs. For additional guidance, see [Technical Note: Covaris g-TUBE DNA Shearing for SMRTbell Prep Kit 3.0 \(102-326-501\)](#) or contact [PacBio Technical Support](#) or your local Field Applications Scientist.

# Example Megaruptor 3 shearing results for a human genomic DNA sample


## Megaruptor 3 System Speed Setting 31



Femto Pulse DNA sizing QC analysis of a human gDNA sample sheared using a Megaruptor 3 with speed setting 31 (1-cycle shear). The fragment size distribution mode is 17.8 kb.

# Example Megaruptor 3 shearing results for microbial genomic DNA samples


## Megaruptor 3 System Speed Setting 40



Femto Pulse DNA sizing QC analyses of four different microbial gDNA samples sheared using a Megaruptor 3 System with speed setting 40 (1-cycle shear). The mean sheared DNA fragment size for all samples is ~7 kb – 10 kb.

## g-TUBE

## 3287 x g [7000 RPM with Eppendorf MiniSpin Plus)



Femto Pulse DNA sizing QC analyses of four different microbial gDNA samples sheared using g-TUBES with a centrifugation speed of 3287 x g. The mean sheared DNA fragment size for all samples is ~7 kb – 10 kb.



# SMRTbell barcoded adapter plate 3.0 is recommended for barcoding WGS samples

For Sequel II and IIe Systems, SMRTbell barcoded adapter plate 3.0 (102-009-200) is available for multiplexing up to 96 microbes per SMRT Cell 8M.

- Contains 96 barcoded adapters to support multiplexed SMRTbell library construction for up to 96 samples using SPK 3.0
- Can be used for Microbial Assembly and any other WGS or amplicon sequencing application that employs barcoded overhang adapters
- Each barcoded adapter contains a 5 bp padding sequence for more uniform ligation performance across different barcode sequences
- Each well on the plate contains a barcoded adapter with a unique 10-base pair PacBio barcode sequence
- Each barcoded adapter is present in only one well and supports a single reaction
- SMRT Link comes pre-installed with the following barcode set FASTA file containing SMRTbell barcoded adapter plate 3.0 barcode sequences\*:  
SMRTbell Barcoded Adapter Plate 3.0 (bc2001-bc2096)

	1	2	3	4	5	6	7	8	9	10	11	12
A	BC 2001	BC 2009	BC 2017	BC 2025	BC 2033	BC 2041	BC 2049	BC 2057	BC 2065	BC 2073	BC 2081	BC 2089
B	BC 2002	BC 2010	BC 2018	BC 2026	BC 2034	BC 2042	BC 2050	BC 2058	BC 2066	BC 2074	BC 2082	BC 2090
C	BC 2003	BC 2011	BC 2019	BC 2027	BC 2035	BC 2043	BC 2051	BC 2059	BC 2067	BC 2075	BC 2083	BC 2091
D	BC 2004	BC 2012	BC 2020	BC 2028	BC 2036	BC 2044	BC 2052	BC 2060	BC 2068	BC 2076	BC 2084	BC 2092
E	BC 2005	BC 2013	BC 2021	BC 2029	BC 2037	BC 2045	BC 2053	BC 2061	BC 2069	BC 2077	BC 2085	BC 2093
F	BC 2006	BC 2014	BC 2022	BC 2030	BC 2038	BC 2046	BC 2054	BC 2062	BC 2070	BC 2078	BC 2086	BC 2094
G	BC 2007	BC 2015	BC 2023	BC 2031	BC 2039	BC 2047	BC 2055	BC 2063	BC 2071	BC 2079	BC 2087	BC 2095
H	BC 2008	BC 2016	BC 2024	BC 2032	BC 2040	BC 2048	BC 2056	BC 2064	BC 2072	BC 2080	BC 2088	BC 2096

Figure illustration of mapping between a specific well location and a unique PacBio barcode sequence on a 96-well plate in the SMRTbell barcoded adapter plate (102-009-200)

Reagent kit quantities support a **single use** of each of the 96 barcoded adapters in the plate for SMRTbell library preparations.

Plate Layout (Excel) [ [Link](#) ]

Barcode Sequences (FASTA) [ [Link](#) ]


Product insert: SMRTbell barcoded adapter plate 3.0 (96 barcodes, 96 samples) [ [Link](#) ]



# Size-selection of SMRTbell libraries with AMPure PB beads is recommended for WGS applications

- AMPure PB Beads are used as the **default size selection method\*** to remove short DNA fragments (<5 kb) and enrich for the long fragments when preparing SMRTbell libraries for whole genome sequencing
- AMPure PB bead size selection of SMRTbell templates is performed follows:
  - Prepare a **35% dilution (v/v)** of the AMPure PB bead stock by adding 1.75 mL of resuspended AMPure PB Beads to 3.25 mL of Elution Buffer (EB). [35% AMPure PB beads solution can be stored at 4°C for 30 days.]
  - Add 3.1X v/v of resuspended, room-temperature 35% AMPure PB beads solution to each sample tube and incubate for 20 min at RT to allow beads to bind to DNA
  - Place sample tubes on a magnetic rack to immobilize AMPure PB beads; wash samples with 80% ethanol 2X; then elute samples in 15  $\mu$ L of EB for 5 min at RT

AMPure PB Beads




With high-quality WGS samples, AMPure PB bead size selection can recover sufficient SMRTbell library material to run **up to ~3 or more SMRT Cells 8M per 3  $\mu$ g of starting input gDNA**


**\* Note:** Although size-selection beads have many advantages, automated DNA size selection systems that utilize gel cassettes offer more flexibility in defining a size cutoff. Three automated DNA size selection tools that may **optionally** be employed for performing size selection on SMRTbell libraries for HiFi WGS applications include the PippinHT, BluePippin and SageELF systems from Sage Science. Note that use of these tools requires **higher input DNA amounts ( $\geq 1.5 \mu\text{g}/\text{SMRT Cell 8M}$ )**. For more information, refer to [Technical Note: Alternative size selection methods for SMRTbell prep kit 3.0 \(TN103-110921\)](#), which provides detailed guidance for size selection of WGS libraries using automated DNA size selection tools or contact [PacBio Technical Support](#) or your local Field Applications Scientist.

# Summary comparison of SPK 3.0 library sample preparation for large genome vs. small genome (microbial / shotgun metagenomic) WGS applications

## Large Genome Workflow



## Small Genome Workflow



Step	Action or metric	Large genome (Single sample)	Small genome (Multiplexed samples)
1	Input DNA quality	<ul style="list-style-type: none"> <li>50% ≥ 30 kb</li> <li>&amp; 90% ≥ 10 kb</li> </ul>	<ul style="list-style-type: none"> <li>90% ≥ 7 kb</li> </ul>
	Input DNA amount*	<ul style="list-style-type: none"> <li>Use ≥1 µg of total input DNA per SMRT Cell 8M</li> </ul>	<ul style="list-style-type: none"> <li>Use ≥300 ng of input DNA per sample, with a total mass ≥1 µg across all samples</li> </ul>
	Megaruptor system shearing	<ul style="list-style-type: none"> <li>Use speed setting 31</li> </ul>	<ul style="list-style-type: none"> <li>Use Speed setting 40</li> </ul>
2	Adapter ligation	<ul style="list-style-type: none"> <li>Use standard (non-barcoded) adapter included with SPK 3.0</li> </ul>	<ul style="list-style-type: none"> <li>Use Barcoded adapter plate 3.0</li> </ul>
3	Size selection	<ul style="list-style-type: none"> <li>Use AMPure PB bead size selection</li> </ul>	<ul style="list-style-type: none"> <li>AMPure PB bead size selection is optional; otherwise perform a standard 1X cleanup with SMRTbell cleanup beads</li> </ul>
4	Pooling	<ul style="list-style-type: none"> <li>N/A</li> </ul>	<ul style="list-style-type: none"> <li>Perform equal mass pooling with barcoded samples</li> </ul>
5	Sequencing prep & data collection	<ul style="list-style-type: none"> <li>30 hr movie time</li> </ul>	<ul style="list-style-type: none"> <li>15 hr movie time</li> </ul>

\* Increase DNA input amounts to ≥1.5 µg per SMRT Cell 8M when using a gel-based automated size selection system.

# How many shotgun metagenomic WGS samples can be multiplexed on a single SMRT Cell 8M?

The overall goals of your project will determine the needed coverage depth

**Question 1:** *What is the estimated abundance of the rarest species you want to observe?*

Example: "I want to see species present at 1% abundance."

→ With 1 SMRT Cell 8M, you can expect ~24,000 HiFi (≥Q20) reads from a 1% abundant species with an 'average' genome size

**Question 2:** *What is your goal?*

In order to achieve...	...You need
Species detection	~100 HiFi reads
Comprehensive gene profiling / discovery*	5-Fold coverage; ~3,000 HiFi reads
Complete genome assembly*	20-Fold coverage; ~12,000 HiFi reads

\* # Reads Needed = Coverage x 5 Mb Genome / 8.5 kb Median HiFi Read Length


# How many shotgun metagenomic WGS samples can be multiplexed on a single SMRT Cell 8M? (cont.)

Example calculation of estimated coverage levels achievable for rare species at different multiplex levels

	1 Sample / SMRT Cell 8M	2 Samples / SMRT Cell 8M	3 Samples / SMRT Cell 8M
Assignable HiFi (≥Q20) Reads per SMRT Cell 8M*	2.4 M	2.4 M	2.4 M
HiFi Reads per Sample	2.4 M	1.2 M	800,000
1% of Reads	24,000 → <b>assembly</b>	12,000 → <b>assembly</b>	8,000 → <b>profiling</b>
0.2% of Reads	4,800 → <b>profiling</b>	2,400 → <b>detection</b>	1,600 → <b>detection</b>

\* Typically, ≥99.5% of HiFi reads have recoverable barcodes.


- The average HiFi read length for metagenomics samples is ~8.5 kb when following our recommended procedure with samples meeting the minimum DNA input quality requirements
- Choose your multiplex level depending on how many reads per rarest-OTU of interest you require for your metagenomic analysis plan



# WGS library sequencing preparation workflow overview

# Sample Setup & Run Design recommendations for WGS libraries

In SMRT Link Sample Setup & Run Design, select 'HiFi Reads', 'Microbial Assembly' or 'Shotgun Metagenomic Profiling or Assembly' for application type



Sequel II binding kit 3.2 & cleanup beads (102-333-300) is recommended for preparing WGS samples for sequencing.

- We recommend using **Sequel II binding kit 3.2 & cleanup beads (102-333-300)** to perform ABC (anneal primer / bind polymerase / clean up complex) with WGS samples
- Refer to **Quick reference card – loading and pre-extension time recommendations for the Sequel II/IIe systems (101-769-100)** for updates to ABC workflow for specific applications

**Sequel II binding kit 3.2 & cleanup beads (102-333-300)** includes the following components:

- Sequencing primer 3.2
- Sequel II polymerase 2.2
- SMRTbell cleanup beads for complex cleanup
- DNA internal control 3.2 (defined 11 kb template bound to Polymerase 2.2)
- Supports  $\geq 24$  binding reactions, and up to 4 SMRT Cells 8M per binding reaction (96 cells total), depending on use case, sample size and concentration

# Sample Setup guidance for WGS samples

Use SMRT Link Sample Setup High-Throughput (HT) mode and follow instructions to perform ABC (anneal primer / bind polymerase / clean up complex) using recommended settings for WGS samples

The screenshot shows the PacBio Sample Setup interface. The 'Version' dropdown is set to 'High-Throughput'. Below the version selector, the text reads 'Sample Setup HT for Sequel II and Sequel IIe'. A table displays one sample setup entry:

<input type="checkbox"/>	Name	Date Created ↓	Number of Samples	Comment	Created By	Locked
<input type="checkbox"/>	Example AAV Sample Setup HT m...	2022-04-12, 01:31:07 PM	2	This batch includes Pooled_AAV_Sample_01 Pooled_AAV_Sample_02	smark	false

The screenshot shows a detailed view of a sample group. The 'Application' field is set to 'HiFi Reads' and is highlighted in green. The 'Binding Kit' field is set to 'Sequel II Binding Kit 3.2' and is also highlighted in green. The 'Cleanup Anticipated Yield' field is set to '75 %' and is highlighted in green. Other parameters include:

- Number of Samples: 2 samples
- SMRT Cells per Sample: 3 cells
- Available Volume per Sample: 15 uL
- Insert Size: 18000 bp
- Sample Concentration: 60 ng/uL
- Recommended Concentration on Plate: 50-90 pM
- Specify Concentration on Plate: 85 pM
- Minimum Pipetting Volume: 1 uL

- **Sample Setup High-Throughput** mode provides a simplified, streamlined workflow to efficiently process either one sample or multiple samples with similar library properties (such as mean insert size and DNA concentration) in parallel
- You can also export the calculated values to a CSV file for **laboratory automation**

**Binding Kit and Cleanup Anticipated Yield fields are auto-filled and highlighted in green after specifying application type**

Example Sample Setup HT mode worksheet for a batch comprised of two human WGS samples.



# Run Design guidance for WGS samples

Follow SMRT Link Run Design instructions to set up a sequencing run using recommended settings for WGS samples

- Choose the appropriate **Application type** for your samples:
  - Select **HiFi Reads** if analyzing large genome WGS samples (default movie collection time = 30 hrs)
  - Select **Microbial Assembly** if analyzing small genome WGS samples (default movie collection time = 15 hrs)
  - Select **Shotgun Metagenomic Profiling or Assembly** if analyzing shotgun metagenomic WGS samples (default movie collection time = 30 hrs)

Include 5mC calls in CpG motifs is enabled by default if **HiFi Reads** or **Custom** is specified as the application type

PacBio Run Design

smark (Lab Tech) ?

Run Information

System Type

SEQUEL II

SEQUEL IIe

Run Name

Example\_Human\_WGS\_Run\_Design

Run Comments

Experiment Name

Experiment ID

Estimated Run Duration (hours): 33.5

Run Reagents / Consumables

1 SMRT Cell

1 sequencing reagent plate

1 mineral oil tube

3 boxes of tips

1 mixing plate

1 sample plate

Sample Information

SAMPLE 1: Human\_WGS\_Library\_01, A01, 30 hour movie, 18000 bp insert

Import from Sample Setup

Select Sample

Application Required

HiFi Reads

Well Sample Name Required

Human\_WGS\_Library\_01

Bio Sample Name Required

HG002

Sample Comment

Sample Well

A01

Template Prep Kit Required

SMRTbell® Prep Kit 3.0

Binding Kit Required

Sequel® II Binding Kit 3.2

Sequencing Kit Required

Sequel® II Sequencing Plate 2.0 (4 rxn)

DNA Control Complex

Sequel® II DNA Internal Control Complex 3.2

Insert Size (bp) Required

18000

Recommended Concentration on Plate (pM)

50-90 pM

On-Plate Loading Concentration (pM) Required

85

Movie Time per SMRT Cell (hours)

30

Use Pre-Extension

YES  NO

Pre-Extension Time (hours)

2

Include 5mC Calls in CpG Motifs

YES  NO

CCS Analysis will be performed on-instrument to produce HiFi .bam files.


Example sample information entered into Run Design for sequencing a human WGS sample.

# Include 5mC calls in CpG motifs option

If selected, kinetic signatures of cytosine bases in CpG motifs will be automatically analyzed to identify the presence of 5mC during on-instrument CCS (Sequel Ii system only) or during CCS analysis in SMRT Link

- **Default setting = YES** when specifying 'HiFi Reads' or 'Custom' application types
- 5mC detection is automatically performed on-instrument with the Sequel Ii system and in SMRT Link with the Sequel II system (data outputs are the **same** for both methods)
- 5mC calls are output in `hifi_reads.bam` as BAM standard MM and ML tags and can be easily visualized in [IGV](#)
- Processing and storage requirements are **minimal**:
  - File size increase is ~5%
  - On-instrument processing time for Sequel Ii systems is ~10 minutes
- Kinetics are not retained in the CCS analysis output by default, but they can **optionally** be retained as well.
- 5mC calls require a **CpG context and symmetrical methylation** (i.e., does not detect hemi-methylated sites)
- Though trained on human data, 5mC detection has been demonstrated to work on non-human data (e.g., plants (Maize)).
- 5mC consensus calling and other tools planned for a future SMRT Link version.
  - For guidance on command line tool options for 5mC analysis, please contact your local PacBio support team or [PacBio Technical Support](#)

Include 5mC Calls in CpG Motifs  YES  NO





**Hypomethylation at transcription start sites**

Example IGV plot demonstrating 5mC detection in HiFi reads for a human HG002 sample. Hypomethylation at active transcription start sites can be easily visualized (unpublished data).

# Include 5mC calls in CpG motifs option (cont.)

- The 5mC CpG detection utility generates the following reports:
  - **CpG Methylation in Reads:** Plots the cumulative percentage of CpG sites in the sample against the predicted probability of methylation. (Report appears in SMRT Link Run QC and Data Management)
  - **CpG Methylation in Reads (Histogram):** Histogram plot displaying the percentage of CpG sites in the sample versus the predicted probability of methylation (Report appears in Data Management only)



# Run Design guidance for multiplexed WGS samples

Run Design setup procedure for automated barcode demultiplexing of WGS library samples barcoded with SMRTbell barcoded adapter plate 3.0

1. Sample is Barcoded: **YES**
2. Barcode Set: Select 'SMRTbell barcoded adapter plate 3.0 (bc2001-bc2096)' (or other custom barcode set if appropriate)
3. Same Barcodes on Both Ends of Sequence: **YES**
4. Assign a **Biological Sample Name** to each barcoded sample using one of two ways: From a (CSV) File or Interactively
5. Specify if barcode demultiplexing is to be performed **on-instrument** (Sequel IIe system only) or in SMRT Link. (Optionally specify Do Not Generate.)

The screenshot shows the PacBio Run Design interface. The 'Sample Information' section is active, and the 'Advanced Options' dropdown is expanded to show 'Barcoded Sample Options'. This section is highlighted with a red box and contains the following configuration:

- 1. Sample Is Barcoded:  YES  NO
- 2. Barcode Set: SMRTbell Barcoded Adapter Plate 3.0
- 3. Same Barcodes on Both Ends of Sequence:  YES  NO
- 4. Assign Bio Sample Names to Barcodes: Interactively
- 5. Demultiplex Barcodes:  ON INSTRUMENT  IN SMRT LINK  DO NOT GENERATE

Buttons for 'Copy', 'Delete', and 'Add Sample' are visible at the bottom of the 'Advanced Options' section.

Example barcoding information entered into Run Design for sequencing a multiplexed WGS sample.



# **WGS data analysis recommendations for variant detection and de novo assembly applications**

# HiFi WGS data analysis recommendations for large genomes

## Using HiFi reads for *de novo* assembly analysis of large genomes

- Perform CCS analysis on-instrument using the Sequel IIe System or in [SMRT Link](#) to generate highly accurate and long single-molecule reads (HiFi reads)
- **10- to 15-fold HiFi read coverage per haplotype** is recommended for most *de novo* assembly projects
  - $Target\ HiFi\ Base\ Yield = [Haploid\ Genome\ Size\ (Gb)] \times [Ploidy\ Level] \times [Target\ HiFi\ Coverage\ per\ Haplotype]$   
E.g., for *de novo* assembly analysis of a 3 Gb diploid genome:  
Recommended Minimum Target HiFi Base Yield = 3 Gb x 2 x 10 = 60 Gb
- Output data in standard file formats, (BAM and FASTA/Q) for seamless integration with downstream analysis tools
- Can use [SMRT Link](#) Genome Assembly analysis application (powered by [IPA](#)) or other third-party software for *de novo* assembly analysis using HiFi reads:
  - [Hifiasm](#)
  - [HiCanu](#)
- Contact PacBio Technical Support ([support@pacb.com](mailto:support@pacb.com)) or your local Bioinformatics Field Applications Scientist for additional information about data analysis recommendations

# HiFi WGS data analysis recommendations for large genomes (cont.)

## Using HiFi reads for variant detection analysis of large genomes

- Perform CCS analysis on-instrument using the Sequel IIe System or in [SMRT Link](#) to generate highly accurate and long single-molecule reads (HiFi reads)
- **≥15-fold HiFi read coverage per sample** is sufficient for most human variant detection projects
  - *Target HiFi Base Yield = [Sample Haploid Genome Size (Gb)] x [Target Coverage per Sample]*  
E.g., For variant detection analysis of a human genome (3 Gb):  
Recommended Minimum Target HiFi Base Yield = 3 Gb x 15 = 45 Gb
- Output data in standard file formats, (BAM and FASTA/Q) for seamless integration with downstream analysis tools
- For detection of small variants (SNVs, InDels <50 bp):
  - Can use third-party software (e.g., Google [DeepVariant](#))
- For detection of structural variants (SVs >50 bp):
  - Can use the Structural Variant Calling application in [SMRT Analysis](#)
- Contact PacBio Technical Support ([support@pacb.com](mailto:support@pacb.com)) or your local Bioinformatics Field Applications Scientist for additional information about data analysis recommendations

# HiFi WGS data analysis recommendations small genomes (microbial multiplexing applications)

## Using HiFi reads for *de novo* assembly and base modification detection analysis of microbial genomes

- Perform CCS analysis on-instrument using the Sequel IIe System or in [SMRT Link](#) to generate highly accurate and long single-molecule reads (HiFi reads)
- **15-fold HiFi read coverage per microbe** is recommended for most *de novo* assembly projects
  - $Target\ HiFi\ Base\ Yield = [Microbe\ Genome\ Size\ (Mb)] \times [Target\ HiFi\ Coverage\ per\ Microbe]$   
E.g., for *de novo* assembly analysis of a 5 Mb microbial genome:  
Recommended Minimum Target HiFi Base Yield = 5 Mb x 15 = 75 Mb
- Output data in standard file formats, (BAM and FASTA/Q) for seamless integration with downstream analysis tools
- Can use [SMRT Link](#) Microbial Genome analysis application for *de novo* assembly and base modification detection analysis using HiFi reads:
  - **Easy to use** (no requirement for laborious parameter input/optimization)
  - **Enables fast and efficient** microbial assembly results using HiFi reads (typical time to result is ~20-60 minutes\* for analysis of a 96-plex microbial data set (up to 375 total sum of genome sizes))
  - **Outputs complete, high-quality** microbial genome assemblies (including chromosomes and plasmids)
- Contact PacBio Technical Support ([support@pacb.com](mailto:support@pacb.com)) or your local Bioinformatics Field Applications Scientist for additional information about data analysis recommendations



# HiFi WGS data analysis recommendations small genomes (microbial multiplexing applications) (cont.)

Use SMRT Link Microbial Genome Analysis application to perform microbial assembly and base modification detection using HiFi reads

PacBio SMRT Analysis

SMRT Analysis / Create New Analysis

1. Select Data 2. Select Analysis

Analysis Application Required

Microbial Genome Analysis

Import Analysis Settings Export

Run Base Modification Analysis

ON  OFF

Find Modified Base Motifs

ON  OFF

Advanced Parameters

- Generate de novo assemblies of small **prokaryotic genomes between 1.9-10 Mb** and companion **plasmids between 2 – 220 kb**, and identify methylated bases and associated nucleotide motifs.
- Optionally include identification of **6mA and 4mC** modified bases and associated **DNA sequence motifs**. (This requires kinetic information.)
  - Unlike 5mC calling, microbial base modification detection is performed **off-instrument** (i.e., in SMRT Link only)
  - This requires a Run Design to specify that **kinetic information be retained** in the CCS analysis output
    - For the Microbial Assembly application type, Run Design **automatically** defaults to specifying **YES** for the '**CCS Analysis Output – Include Kinetics Information**' field

Application Required


Microbial Assembly

CCS Analysis Output – Include Kinetics Information  YES  NO

**Note:** This combines and **replaces** the Microbial Assembly and Base Modification Analysis applications in SMRT Link releases prior to v11.0.

# HiFi WGS data analysis recommendations small genomes (microbial multiplexing applications) (cont.)

View SMRT Link microbial assembly results, detected base modifications and identified modified base motifs in a single analysis job report



# HiFi WGS data analysis recommendations small genomes (shotgun metagenomic applications)

## HiFi reads are compatible with third-party metagenomics data analysis tools

- Perform CCS analysis on-instrument using the Sequel IIe System or in [SMRT Link](#) to generate highly accurate and long single-molecule reads (HiFi reads)
- Output data in standard file formats, (BAM and FASTA/Q) for seamless integration with downstream analysis tools
- Use [PacBio metagenomics tools](#) available on GitHub for **taxonomic classification and functional gene profiling** using HiFi reads
- Can perform **metagenomic shotgun assembly** directly with HiFi reads using [Hifiasm](#) and evaluate & extract **metagenome-assembled genomes** using PacBio [HiFi-MAG-Pipeline](#) tool available on GitHub (see Portik *et al.*)
- Contact PacBio Technical Support ([support@pacb.com](mailto:support@pacb.com)) or your local Bioinformatics Field Applications Scientist for additional information about data analysis recommendations




# WGS library example sequencing performance data


# Example library QC and sequencing performance for human WGS libraries prepared with SMRTbell prep kit 3.0

## SMRTbell library QC and primary sequencing metrics

### Size-Selected Library QC




### Raw Data Report



A secondary left hand peak may also be visible depending on DNA sample quality

### CCS Analysis Report



Input gDNA for Megaruptor 3 shearing	3000 ng
Post-shearing recovery (%)*	2620 ng (87%)
Final yield of AMPure PB bead Size-selected library (%)**	1070 ng (36%)

\* Post-shearing recoveries typically ranged from ~70% to >95% when using input human DNA samples (1 µg to 5 µg)

\*\* Final post-size selected library yields typically ranged from ~25% to ~50% when using input human DNA samples (1 µg to 5 µg)

Raw Base Yield	617.65 Gb
Mean Polymerase Read Length	102.8 kb
P0	23.4%
P1	75.0%
P2	1.6%

Example sequencing metrics for a human WGS sample run with Binding Kit 3.2 (Polymerase 2.2) / 85 pM on-plate concentration / 30-h movie time / 2-h Pre-Extension Time / Adaptive Loading Target = 0.85


HiFi Reads	2.7 M
HiFi Base Yield	39.2 Gb
Mean HiFi Read Length	14,490 bp
Median HiFi Read Quality	Q34
HiFi Read Mean # of Passes	12

For SPK 3.0 human WGS libraries, per-SMRT Cell HiFi base yields typically ranged from ~28 to 39 Gb.


# Example sequencing performance for human WGS libraries prepared with SMRTbell prep kit 3.0 (cont.)

## Primary sequencing metrics (Cont.)

### Yield by HiFi Read Length



### Read Quality Distribution



### HiFi Read Length Summary

Read Length (bp)	Reads	Reads (%)	Yield (bp)	Yield (%)
≥ 0	2,707,732	100	39,236,168,651	100
≥ 5,000	2,664,322	98	39,051,919,399	100
≥ 10,000	2,353,137	87	36,541,368,326	93
≥ 15,000	1,164,272	43	21,435,305,025	55
≥ 20,000	294,460	11	6,522,779,501	17
≥ 25,000	21,062	1	559,040,421	1
≥ 30,000	1,012	0	35,294,569	0
≥ 35,000	388	0	15,240,023	0
≥ 40,000	129	0	5,578,841	0

### HiFi Read Quality Summary

Read Quality (Phred)	Reads	Reads (%)	Yield (bp)	Yield (%)
≥ Q20	2,707,732	100	39,236,168,651	100
≥ Q30	1,811,377	67	25,413,473,886	65
≥ Q40	679,582	25	8,150,599,400	21
≥ Q50	146,257	5	1,355,549,531	3

# Example *de novo* assembly performance for human WGS libraries prepared with SMRTbell prep kit 3.0

HiFi WGS data sets generated with SPK 3.0 provide highly contiguous and highly accurate assemblies

HG002 Library ID	Contig_N50_Mbp
64009e_s10_cov30	35.4
64012e_s10_cov30	36.5
64015e_s10_cov30	36.7
64438e_s10_cov30	34.1
64441e_s10_cov30	33.6

HG002 Library ID	deNovo_asm_QV
64009e_s10_cov30	48.3
64012e_s10_cov30	48.3
64015e_s10_cov30	48.2
64438e_s10_cov30	48.3
64441e_s10_cov30	48.2

- Data were generated from five different human HG002 WGS libraries run on five different Sequel IIe system instruments
- Data were subsampled to 30-fold coverage and assembled using SMRT Link Genome Assembly analysis application

# Example variant detection performance for human WGS libraries prepared with SMRTbell prep kit 3.0

HiFi WGS data sets generated with SPK 3.0 provide highly accurate variant calls

HG002 Library ID	INDEL.F1_Score
64009e_s10_cov30	0.995
64012e_s10_cov30	0.994
64015e_s10_cov30	0.994
64438e_s10_cov30	0.993
64441e_s10_cov30	0.994

HG002 Library ID	SNP.F1_Score
64009e_s10_cov30	0.999
64012e_s10_cov30	0.999
64015e_s10_cov30	0.999
64438e_s10_cov30	0.999
64441e_s10_cov30	0.999

- Data were generated from five different human HG002 WGS libraries run on five different Sequel IIe system instruments
- Data were subsampled to 30-fold coverage and analyzed with DeepVariant



# Example sequencing performance for a 96-plex microbial WGS library prepared with SMRTbell prep kit 3.0

## Sample preparation workflow

### Experiment design

- 24 different microbes; each ligated independently to 4 different barcodes for 96-plex
- Selected microbes relevant to food safety and human health represent a range of genome sizes, GC content, and plasmid composition
- Total sum of genome sizes = 375 Mb

### SMRTbell library construction

- 1 µg of input gDNA per microbe for shearing
- Target shear size = 7 kb – 10 kb
- 500 ng of sheared DNA into library prep
- Symmetrically barcoded samples using SMRTbell Barcoded Adapter Plate 3.0 Kit (102-009-200)
- No size-selection performed

#### Microbial genome assembly complexity\*

**Class I** – Have few repeats except for the rDNA operon sized 5 to 7 kb

**Class II** - Class II genomes have many repeats, such as insertion sequence elements, but none greater than 7 kb.

**Class III** - Contain large, often phage-related, repeats >7 kb.


Microbial species	Genome size (bp)	GC content (%)	Microbial genome complexity	Barcode names
<i>Acinetobacter baumannii</i> AYE	3,960,239	39.35	Class 3	bc2001 / bc2025 / bc2049 / bc2073
<i>Bacillus cereus</i> 971	5,430,163	35.29	Class 1	bc2002 / bc2026 / bc2050 / bc2074
<i>Bacillus subtilis</i> W23	4,045,592	43.5	Class 1	bc2003 / bc2027 / bc2051 / bc2075
<i>Burkholderia cepacia</i> UCB 717	8,569,621	66.6	Class 3	bc2004 / bc2028 / bc2052 / bc2076
<i>Burkholderia multivorans</i> 249	7,008,277	66.68	Class 3	bc2005 / bc2029 / bc2053 / bc2077
<i>Enterococcus faecalis</i> OG1RF	2,739,503	37.75	Class 1	bc2006 / bc2030 / bc2054 / bc2078
<i>Escherichia coli</i> H10407	5,393,109	50.71	Class 1	bc2007 / bc2031 / bc2055 / bc2079
<i>Escherichia coli</i> K12 MG1655	4,642,522	50.79	Class 1	bc2008 / bc2032 / bc2056 / bc2080
<i>Helicobacter pylori</i> J99	1,645,141	39.19	Class 1	bc2009 / bc2033 / bc2057 / bc2081
<i>Klebsiella pneumoniae</i> BAA-2146	5,780,684	56.97	Class 2	bc2010 / bc2034 / bc2058 / bc2082
<i>Listeria monocytogenes</i> Li2	2,950,984	37.99	Class 1	bc2011 / bc2035 / bc2059 / bc2083
<i>Listeria monocytogenes</i> Li23	2,979,685	38.19	Class 1	bc2012 / bc2036 / bc2060 / bc2084
<i>Methanocorpusculum labreanum</i> Z	1,804,962	50.5	Class 1	bc2013 / bc2037 / bc2061 / bc2085
<i>Neisseria meningitidis</i> FAM18	2,194,814	51.62	Class 3	bc2014 / bc2038 / bc2062 / bc2086
<i>Neisseria meningitidis</i> Serogroup B	2,304,579	51.44	Class 1	bc2015 / bc2039 / bc2063 / bc2087
<i>Rhodopseudomonas palustris</i> CGA009	5,459,213	64.9	Class 3	bc2016 / bc2040 / bc2064 / bc2088
<i>Salmonella enterica</i> LT2	4,950,860	52.24	Class 1	bc2017 / bc2041 / bc2065 / bc2089
<i>Salmonella enterica</i> Ty2	4,791,947	52.05	Class 1	bc2018 / bc2042 / bc2066 / bc2090
<i>Staphylococcus aureus</i> Seattle 1945	2,806,348	32.86	—	bc2019 / bc2043 / bc2067 / bc2091
<i>Staphylococcus aureus</i> USA300_TCH1516	2,872,915	32.7	Class 1	bc2020 / bc2044 / bc2068 / bc2092
<i>Streptococcus pyogenes</i> Bruno	1,844,942	38.48	—	bc2021 / bc2045 / bc2069 / bc2093
<i>Thermanaerovibrio acidaminovorans</i> DSM6589	1,852,980	63.78	Class 1	bc2022 / bc2046 / bc2070 / bc2094
<i>Treponema denticola</i> A	2,842,721	37.87	—	bc2023 / bc2047 / bc2071 / bc2095
<i>Vibrio parahaemolyticus</i> EB101	5,146,979	45.33	Class 1	bc2024 / bc2048 / bc2072 / bc2096

\* Koren, S. et al. (2013) Reducing assembly complexity of microbial genomes with single-molecule sequencing. [Genome Biol 14, R101](#)


# Example sequencing performance for A 96-plex microbial WGS library prepared with SMRTbell prep kit 3.0 (cont.)

## Primary sequencing metrics


### Size-Selected Library QC



### Raw Data Report



### CCS Analysis Report



Input gDNA per microbe For Megaruptor 3 shearing	1 µg
Input sheared DNA per microbe For library construction	500 ng
Mean SMRTbell library construction yield per microbe before pooling (%)*	32%

\* LC yields ranged from 17 – 52% across 96 microbes. Microbial libraries were barcoded with SMRTbell Barcoded Adapter Plate 3.0 Kit, independently purified with SMRTbell Cleanup Beads after nuclease treatment, and then pooled for sequencing on a single SMRT Cell 8M.

Raw Base Yield	224.94 Gb
Mean Polymerase Read Length	79.2 kb
P0	63.3
P1	35.6
P2	1.1


90 pM on-plate concentration / 15-h movie time / 2-h Pre-Extension Time / Adaptive Loading Target = 0.85

HiFi Reads	1.8 M
HiFi Base Yield	13.8 Gb
Mean HiFi Read Length	7,881 bp
Median HiFi Read Quality	Q38
HiFi Read Mean # of Passes	14

# Example sequencing performance for A 96-plex microbial WGS library prepared with SMRTbell prep kit 3.0 (cont.)

## Primary sequencing metrics (Cont.)


### HiFi Read Length Distribution



### HiFi Read Length Summary

Read Length (bp)	Reads	Reads (%)	Yield (bp)	Yield (%)
≥ 0	1,756,336	100	13,842,769,905	100
≥ 5,000	1,451,205	83	12,852,531,184	93
≥ 10,000	416,073	24	4,979,117,856	36
≥ 15,000	24,686	1	397,536,526	3
≥ 20,000	206	0	4,478,104	0
≥ 25,000	13	0	350,714	0
≥ 30,000	2	0	61,197	0
≥ 35,000	0	0	0	0
≥ 40,000	0	0	0	0

### Read Quality Distribution





### HiFi Read Quality Summary

Read Quality (Phred)	Reads	Reads (%)	Yield (bp)	Yield (%)
≥ Q20	1,756,336	100	13,842,769,905	100
≥ Q30	1,411,586	80	10,691,739,197	77
≥ Q40	687,603	39	4,137,789,984	30
≥ Q50	302,228	17	1,376,190,958	10

# Example sequencing performance for A 96-plex microbial WGS library prepared with SMRTbell prep kit 3.0 (cont.)

## Barcode demultiplexing results


Value	Analysis Metric
96	Unique Barcodes
1,731,704	Barcoded Reads
18,038	Mean Reads
64,709	Max. Reads
4,565	Min. Reads
7,856	Mean Read Length
24,632	Unbarcoded Reads
98.66%	Percent Bases in Barcoded Reads
98.59%	Percent Barcoded Reads



- All 96 barcodes detected
- Mean # of barcoded HiFi reads per microbe is ~18,000
- Mean HiFi base coverage per microbe is 36-fold (Range is 19- to 63-fold)

# Example sequencing performance for A 96-plex microbial WGS library prepared with SMRTbell prep kit 3.0 (cont.)

## HiFi *de novo* assembly results – assembled chromosomes




- Achieved 1 Contig / Chromosome for 92 out of 96 assemblies
- For all 96 microbes, most or all chromosomal assemblies were complete and of the expected sizes

Microbial assembly statistics from a 96-plex pool of bacteria relevant to food safety and human health. These data were generated on the Sequel II system and assembled with the fully automated HiFi-based Microbial Assembly application in SMRT Link using the default parameters, without any manual curation. [Download](#) and explore the data yourself.

# Example sequencing performance for A 96-plex microbial WGS library prepared with SMRTbell prep kit 3.0 (cont.)

HiFi *de novo* assembly results – representative assembly accuracies



- Accuracy of representative samples from a 96-plex microbial whole genome sequencing run on a Sequel II System
- With HiFi data and the Microbial Assembly application in SMRT Link, **genome assemblies are consistently >99.99% accurate**



# Technical documentation & applications support resources

# Technical resources for WGS library preparation, sequencing & data analysis

## DNA extraction literature

- Circulomics Nanobind [Application notes](#)
- Circulomics Nanobind [Kit handbooks](#)
- Circulomics Nanobind [Protocols](#)
- Technical note: Sample preparation for PacBio HiFi sequencing from human whole blood ([102-326-500](#))
- Technical note: Preparing DNA for PacBio HiFi sequencing – Extraction and quality control ([102-193-651](#))
- Technical note: Preparing samples for PacBio whole genome sequencing for *de novo* assembly – collection and storage ([TN100-040518](#))

## Sample preparation literature

- Application brief: Whole genome sequencing (WGS) for *de novo* assembly – Best practices ([102-193-627](#))
- Application brief: Variant detection using whole genome sequencing with HiFi reads – Best practices ([102-193-604](#))
- Overview – Sequel systems application options and sequencing recommendations ([101-851-300](#))
- Procedure & checklist – Preparing whole genome and metagenome sequencing libraries using SMRTbell prep kit 3.0 ([102-166-600](#))
- Quick reference card – Loading and pre-extension recommendations for the Sequel II and IIe systems ([101-769-100](#))
- Technical note: Alternative size selection methods for SMRTbell prep kit 3.0 ([TN103-110921](#))
- Technical note: Covaris g-TUBE DNA shearing for SMRTbell prep kit 3.0 ([102-326-501](#))
- Technical overview – Whole genome and metagenome library preparation using SMRTbell prep kit 3.0 ([102-390-900](#))



# Technical resources for WGS library preparation, sequencing & data analysis (cont.)

## Data analysis resources

- SMRT Link v11.0 user guide ([102-278-200](#))
- SMRT Tools v11.0 reference guide ([102-278-500](#))
- Sequel II and IIe systems: Data files ([102-144-100](#))

## Example PacBio data sets

Whole genome sequencing application	Dataset	Data type	PacBio system
Assembly and variant detection	<a href="#">Homo sapiens – GIAB sample HG002</a>	HiFi Reads	Sequel II System
Assembly	<a href="#">Oryza sativa – MH63</a>	HiFi Reads	Sequel IIe System
Assembly (low-DNA input)	<a href="#">Anopheles gambiae – 2 plex</a>	HiFi Reads	Sequel II System
Assembly (ultra-low DNA input)	<a href="#">Phlebotomus papatasi, Homo sapiens, Drosophila melanogaster</a>	HiFi Reads	Sequel II System
Assembly	<a href="#">Food safety and infectious microbes – 96 plex</a>	HiFi Reads	Sequel II System
5mC detection at CpG sites	<a href="#">Human HG002 CpG methylation status</a>	HiFi Reads	Sequel IIe System

# Technical resources for WGS library preparation, sequencing & data analysis (cont.)

## Publications

- Nurk S. et al. (2022) The complete sequence of a human genome. Science. 376:44-53 doi: [10.1126/science.abj6987](https://doi.org/10.1126/science.abj6987)
- Noyes, D.N. et al. (2022) Familial long-read sequencing increases yield of de novo mutations. American journal of human genetics. 109:631-646. doi: <https://doi.org/10.1016/j.ajhg.2022.02.014>
- Bickhart, D.M. et al. (2022) Generating lineage-resolved, complete metagenome-assembled genomes from complex microbial communities. Nature biotechnology. doi: [10.1038/s41587-021-01130-z](https://doi.org/10.1038/s41587-021-01130-z)
- Lefoulon, E. et al. (2021) Greenhead (Tabanus nigrovittatus) Wolbachia and its microbiome: A preliminary study. Microbiol Spectr. 9(2):e0051721. doi: [10.1128/Spectrum.00517-21](https://doi.org/10.1128/Spectrum.00517-21)

## Webinars

- PacBio Webinar (2022): Unlocking the genome with long-read sequencing in genetic disease research [ [Link](#) ]
- PacBio Webinar (2021): Methylation detection with PacBio HiFi sequencing [ [Link](#) ]
- PacBio Webinar (2021): HiFi sequencing: see what you've been missing [ [Link](#) ]
- PacBio Webinar (2021): Integrated rare disease using long-read genome sequencing [ [Link](#) ]
- PacBio Webinar (2021): Getting the most out of your breeding program with DNA and RNA sequencing [ [Link](#) ]
- PacBio Webinar (2021): Generation of lineage-resolved complete metagenome-assembled genomes in complex microbial communities [ [Link](#) ]



# **APPENDIX: Genomic DNA isolation & QC recommendations for PacBio WGS library sample preparation**

Sequel IIe System ICS v11.0 / SMRT Link v11.0

# Genomic DNA isolation & QC recommendations for PacBio WGS sample preparation

## Technical contents

1. Sample collection, preparation, and storage for SMRT sequencing
2. Genomic DNA extraction, QC and handling for SMRT sequencing
3. Example protocols and kit solutions for high-molecular weight genomic DNA isolation
4. Methods for evaluation of genomic DNA quality
5. Cleanup of genomic DNA and SMRTbell libraries
6. Storage and shipping of genomic DNA and SMRTbell libraries
7. DNA sample extraction literature resources



# Sample collection, preparation, and storage for SMRT sequencing

# Sample collection, preparation, and storage for SMRT sequencing whole genome sequencing projects

To obtain the highest quality genomic DNA, it is important to start with sample types compatible with high molecular weight (HMW) DNA extraction methods

## PacBio Technical note: Sample prep ([TN100-040518](#))

- Provides general guidance on biological sample collection, preparation, and storage across a range of commonly encountered sample types including Vertebrates, Invertebrates, Arthropods, Fungi and Plants

### Invertebrates

When sampling from invertebrates, it is recommended to use a single individual. Some invertebrates have mucous membranes that inhibit the ability to obtain high-quality DNA. Please consider an extra cleanup step of the DNA if mucous-coated samples are used (see cleanup protocol [here](#)). It is also common to encounter invertebrate samples that are not easily separated from contaminants or do not have cell-dense tissues readily available. In these cases, sperm can be used as the input sample.

	Sample Type	Sample Storage
1	Cell-dense tissue (brain, kidney, muscle, etc.)	Fresh (within ~24 hours of collection and kept cold) or flash frozen with liquid nitrogen and stored at -80°C should be viable for several months.
2	Cell-culture	Room temperature cell-suspension or pellet is best entering the DNA isolation step. However, cryopreserved cells (slow frozen in cryoprotectant) should be viable at -80°C for several months.
3	Sperm	Fresh collected, room temperature sperm samples are best entering the DNA isolation step. However, customers have been successful with sperm frozen at -20°C.

**Table 2 - Invertebrate sample types in order of preference**

**Note:** It is NOT recommended to use liver tissue as input sample due to high abundance of enzymes that may degrade DNA.



PacBio Technical note [TN100-040518](#): Preparing samples for PacBio whole genome sequencing for *de novo* assembly – Collection and storage

## Circulomics Nanobind kit handbooks

[Nanobind Kit Handbooks](#) contain information on HMW DNA isolation kit specifications, general tips, tissue preservation recommendations, extraction and sequencing performance, and troubleshooting tips.

- Nanobind CBB Big DNA Kit Handbook ([HBK-CBB-001](#))
- Nanobind Plant Nuclei Big DNA Kit Handbook ([HBK-PLT-001](#))
- Nanobind Tissue Big DNA Kit Handbook ([HBK-TIS-001](#))

### RNAIater-Preserved Tissues

Tissues that are preserved in RNAIater prior to freezing or storage should have excess RNAIater solution removed. After placing the tissue on a clean, chilled surface, wick away excess RNAIater liquid using a Kimwipe.

### Ethanol-Preserved Tissues

Tissues that are preserved in ethanol prior to freezing or storage require pre-treatment before extraction to remove the ethanol.

#### 1. Prepare EtOH Removal Buffer

- 400 mM NaCl

Nanobind Tissue Big DNA Kit Handbook

Document ID: HBK-TIS-001

Release Date: 03/24/2021

For extraction of HMW (50 kb – 300+ kb) and UHMW (50 kb – 1+ Mb) genomic DNA from tissue



Visit the Circulomics [Nanobind Support](#) website to find the latest resources for using Nanobind Kits for HMW DNA extraction and recommended tissue preservation methods.



# Genomic DNA extraction, QC and handling for SMRT sequencing




# DNA extraction, QC and handling for SMRT sequencing whole genome sequencing projects

Starting with high-quality, high molecular weight (HMW) genomic DNA (gDNA) will result in longer libraries and better performance during sequencing

## PacBio Technical note: DNA prep ([102-193-651](#))

- Provides recommendations, tips and tricks for the extraction of genomic DNA, as well as assessing and preserving the quality and size of your DNA sample to be used for PacBio HiFi sequencing for *de novo* assembly
- Topics covered include:
  - DNA extraction
    - Commercially available kits across a wide variety of input sample types
    - Resource for alternative DNA extraction methods
  - DNA quality control (QC)
    - DNA quantification, purity, size, and damage
    - Use of nucleic acid stabilizers
    - DNA storage and shipping
  - General best practices for DNA extraction for PacBio sequencing
- This technical note also includes an example dataset for a California Redwood tree DNA sample that was isolated using the **Nanobind Plant Nuclei Big DNA Kit** ([NB-900-801-01](#)) from [Circulomics](#)



PacBio Technical Note: Preparing DNA for PacBio HiFi sequencing - Extraction and quality control ([102-193-651](#))

# DNA extraction, QC and handling for SMRT sequencing whole genome sequencing projects (cont.)

Whole blood is a common and easily accessible source of DNA that – with proper handling – provides high-quality input for PacBio HiFi sequencing

## PacBio Technical Note: Sample preparation for PacBio HiFi sequencing from human whole blood ([102-326-500](#))

- To define the best practices for handling human whole blood samples, we tested the effect of anticoagulant, sample storage time, storage conditions, and white blood cell count on the sequencing performance of DNA extracted using the **Nanobind CBB Big DNA Kit** ([NB-900-001-01](#)).

Stage	Variable	Best practice for PacBio HiFi sequencing
Before DNA extraction	Sample type	Human whole blood
	Anticoagulant	Potassium EDTA (K <sub>2</sub> EDTA)
	Sample storage temperature	4 ± 3°C
	Sample storage time	≤ 2 days from collection to extraction
DNA extraction	Volume of whole blood	200 µL
	White blood cell (WBC) count	≥ 4 × 10 <sup>6</sup> cells/mL for ≥ 3 µg of DNA
	DNA extraction method	Nanobind CBB Big DNA kit
After DNA extraction	DNA storage	Rest 1 day at ambient temperature, then store at 4 ± 3°C
	DNA size distribution	<ul style="list-style-type: none"> <li>90% of DNA ≥ 10 kb (genomic quality number at 10 kb ≥ 9.0)</li> <li>50% of DNA ≥ 30 kb (genomic quality number at 30 kb ≥ 5.0)</li> </ul>
	UV absorbance	<ul style="list-style-type: none"> <li>A260/280 nm ≥ 1.7</li> <li>A260/230 nm ≥ 1.5</li> </ul>

Technical note  
**SAMPLE PREPARATION FOR PACBIO HIFI SEQUENCING FROM HUMAN WHOLE BLOOD**

**Introduction**  
 Whole blood is a common and easily accessible source of DNA that – with proper handling – provides high-quality input for PacBio HiFi sequencing. To define the best practices for handling human whole blood samples, we tested the effect of anticoagulant, sample storage time, storage conditions, and white blood cell count on the sequencing performance of DNA extracted using the Nanobind® CBB Big DNA kit.

**Summary**

Stage	Variable	Best practice for PacBio HiFi sequencing
Before DNA extraction	Sample type	Human whole blood
	Anticoagulant	Potassium EDTA (K <sub>2</sub> EDTA)
	Sample storage temperature	4 ± 3°C
	Sample storage time	≤ 2 days from collection to extraction
DNA extraction	Volume of whole blood	200 µL
	White blood cell (WBC) count	≥ 4 × 10 <sup>6</sup> cells/mL for ≥ 3 µg of DNA
	DNA extraction method	Nanobind CBB Big DNA kit
After DNA extraction	DNA storage	Rest 1 day at ambient temperature, then store at 4 ± 3°C
	DNA size distribution	<ul style="list-style-type: none"> <li>90% of DNA ≥ 10 kb (genomic quality number at 10 kb ≥ 9.0)</li> <li>50% of DNA ≥ 30 kb (genomic quality number at 30 kb ≥ 5.0)</li> </ul>
	UV absorbance	<ul style="list-style-type: none"> <li>A260/280 nm ≥ 1.7</li> <li>A260/230 nm ≥ 1.5</li> </ul>

**PacBio**

PacBio Technical Note: Sample preparation for PacBio HiFi sequencing from human whole blood ([102-326-500](#))

# DNA extraction, QC and handling for SMRT sequencing whole genome sequencing projects (cont.)

[Circulomics](#) offers innovative products for HMW DNA extraction to sequence common samples as well as more challenging samples such as animal tissue, insects, fungi, and plants

## Circulomics Nanobind kit handbooks

- [Nanobind Kit Handbooks](#) contain information on HMW DNA isolation kit specifications, general tips, tissue preservation recommendations, extraction and sequencing performance, and troubleshooting tips.

## Circulomics Nanobind protocols

- [Nanobind Protocols](#) posted on the Circulomics website are always our most up to date versions and should take precedence over the Kit Handbooks. If an Application Note exists, that should be consulted first.

## Circulomics Nanobind application notes

- [Nanobind Application Notes](#) provide detailed protocols for challenging or interesting samples with supporting sequencing data.

## Circulomics Nanobind developmental protocols

- Developmental protocols exist for many sample types such as insect, worms, and fungal samples (contact [Circulomics](#) for details)

Nanobind CBB Big DNA Kit Handbook

Document ID: HBK-CBB-001

Release Date: 3/24/2021

For extraction of HMW (50 kb – 300+ kb) and UHMW (50 kb – 1+ Mb) genomic DNA from cultured cells, cultured bacteria, and blood



Visit the Circulomics [Nanobind Support](#) website to find the latest resources for using Nanobind Kits for HMW DNA extraction.

# General recommendations for isolating high-molecular weight (HMW) genomic DNA

## Before gDNA extraction:

- Use fresh or flash-frozen tissue
- Store flash-frozen tissue at -80°C and avoid freeze-thaw cycles
- Do not store blood samples longer than a few days at 4-8°C before DNA extraction
- Microbial gDNA Isolation:
  - ❑ Avoid culture incubation in complex or rich media
  - ❑ Harvesting from several replicate cultures rather than a single, high-density culture is preferred
  - ❑ Extraction of small culture volumes is preferred over large volumes to avoid accumulating high concentrations of potentially inhibiting secondary components

## During gDNA extraction:

- Mechanically disrupt tissues using TissueRuptor (QIAGEN), Dounce homogenizer or liquid nitrogen grinding
- Inactivate nucleases and DNA binding proteins with a protease, such as proteinase K
- Remove all RNA with RNase A
- It is preferable to avoid oxidative agents such as phenol and/or chloroform if possible to minimize DNA damage
- Resuspend, or elute, DNA in a low salt buffer, such as 10 mM Tris-HCl pH 8.0-9.0 + 0.1 mM EDTA

# General recommendations for isolating high-molecular weight (HMW) genomic DNA (cont.)

## After gDNA extraction:

- Check concentration on both the NanoDrop and Qubit systems for concordance
- High-quality, pure DNA typically shows a A260/280 ratio  $\geq 1.8$  and a A260/230 ratio  $\geq 2.0$
- If gel purification is required, avoid using ethidium/UV based visualization methods. One alternative is to use SYBR Safe (Invitrogen) and visualize with blue light
- To help resuspend HWM DNA, pipette mix 1-10 times with a standard P200 pipette tip. Allow the DNA to rest overnight at 25°C.
- Overheating can introduce DNA damage. Inactivate DNase as recommended by the vendor kit. It is best to avoid heat inactivation when possible
- DNA storage conditions: 4°C (short-term); -20°C / -80°C (long-term)
- Ideally proceed to SMRTbell library preparation with freshly isolated DNA whenever possible





# Example protocols and kit solutions for high-molecular weight genomic DNA isolation

# Example high-molecular weight genomic DNA isolation protocols and kit solutions

Cultured cells, cultured bacteria (gram negative / gram positive), and blood

## Recommended kit product or protocol

### **Circulomics Nanobind CBB Big DNA Kit** ([NB-900-001-01](#))

- The Nanobind CBB Big DNA Kit is designed for rapid extraction of really big DNA from cultured cells, cultured bacteria (gram negative and gram positive), and blood (nucleated and non-nucleated).
- Includes protocols for HMW DNA (50 kb - 300+ kb) and UHMW DNA (50 kb - 1 Mb+) extraction.
- Typical yields are 5 µg - 100+ µg per extraction depending on input.
- A high-throughput version of this kit is available for automated extractions on Thermo Fisher KingFisher instruments. (Contact [Circulomics](#) for details.)



## Other kit products or protocols\*

QIAGEN Genomic-tip 20/100/500/G Kit [ [Link](#) ]

QIAGEN Genra Puregene Kit [ [Link](#) ]

Geneaid DNA Isolation (Bacteria) Kit [ [Link](#) ]

Lucigen Masterpure Kit [ [Link](#) ]

JGI Bacterial Genomic DNA Isolation Protocol Using CTAB [ [Link](#) ]

NEB Monarch HMW DNA Extraction Kit [ [Link](#) ]

## Plant Tissue

### Recommended kit product or protocol

#### ***Circulomics Nanobind Plant Nuclei Big DNA Kit*** ([NB-900-801-01](#))

- The Nanobind Plant Nuclei Big DNA Kit is designed for rapid extraction of really big DNA from plant nuclei.
- First, nuclei are isolated from 1–5 g of plant tissue using one of the recommended nuclei isolation protocols. Then, HMW DNA is extracted from the nuclei using Nanobind disks. Each of the two purification steps (i.e., nuclei isolation + Nanobind extraction) removes different contaminants from the sample, resulting in clean, HMW DNA from even the most challenging plant species.
- For many plant types, big DNA (up to 300+ kb) can be obtained in <1 hour from the nuclei pellet stage.
- Typical yields are 5 µg - 20+ µg per extraction depending on input



### Other kit products or protocols\*

*QIAGEN Genomic-tip 20/100/500/G Kit* [ [Link](#) ]

*Unsupported Protocol – Switchgrass (*Panicum virgatum*) DNA isolation [USDA]* [ [Link](#) ]

*Unsupported Protocol – DNA extraction of *Chlamydomonas* using CTAB [JGI]* [ [Link](#) ]

*QIAGEN User-Developed Protocol: Isolation of genomic DNA from plants and filamentous fungi using the QIAGEN Genomic-tip Kit* [ [Link](#) ]

*Modified QIAGEN Genomic-tip Protocol [King Abdullah University of Science and Technology]* [ [Link](#) ]



## Animal Tissue

### Recommended kit product or protocol

#### ***Circulomics Nanobind Tissue Big DNA Kit (NB-900-801-01)***

- The Nanobind Tissue Big DNA Kit is designed for rapid extraction of really big DNA from diverse animal tissues including mammalian, fish, avian, mollusk, and crustacean samples.
- Includes protocols for HMW DNA (50 kb - 300+ kb) and UHMW DNA (50 kb - 1+ Mb) extraction.
- Typical yields are 5 ug - 100+ ug per extraction depending on tissue type and input.
- Note: Insect, worm and fungal samples are not officially supported at this time. However, select insects, worms and fungi can be processed with supplemental buffers. Contact [Circulomics](#) for more details.



### Other kit products or protocols\*

QIAGEN Genomic-tip 20/100/500/G Kit [ [Link](#) ]

QIAGEN Genra Puregene Kit [ [Link](#) ]

QIAGEN User-Developed Prototocol: Isolation of genomic DNA from mosquitoes or other insects using the QIAGEN Genomic-tip Kit [ [Link](#) ]

QIAGEN User-Developed Prototocol: Purification of archive-quality DNA from 10–20 mg fish tissue using the Genra Puregene Tissue Kit or Genra Puregene Mouse Tail Kit [ [Link](#) ]

Unsupported Protocol – Genra Puregene Cell Kit (Qiagen) DNA Isolation [Univ. Washington][ [Link](#) ]

Macherey-Nagel™ NucleoBond™ AXG 20/100/500 Gravity-flow Columns [ [Link](#) ]

## Yeast and Fungi

### Recommended kit product or protocol

- Contact [Circulomics](#) for details on how to process select yeast and fungal samples with Nanobind Kits using supplemental buffers.



**Technical Support:**  
[support@circulomics.com](mailto:support@circulomics.com)



### Other kit products or protocols\*

QIAGEN Genomic-tip 20/100/500/G Kit [ [Link](#) ]

QIAGEN Genra Puregene Kit [ [Link](#) ]

QIAGEN DNeasy PowerLyzer PowerSoil Kit [ [Link](#) ]

GeneJET Plant Genomic DNA Purification Kit [ [Link](#) ]

Zymo Research Fungal/Bacterial DNA MidiPrep Kit [ [Link](#) ]

JGI Yeast and Fungal DNA Isolation Protocol [ [Link](#) ]

*Have challenging or interesting samples to process?*

Please contact [Circulomics](https://www.circulomics.com) for recommendations on how to process your challenging or interesting samples using Nanobind Kits.



**Technical Support:**  
[support@circulomics.com](mailto:support@circulomics.com)






# Cleanup of genomic DNA and SMRTbell libraries

# General recommendations for cleanup of genomic DNA and SMRTbell libraries

## AMPure PB beads can be used for cleanup of genomic DNA

- A quick and very effective check for sample purity is to compare the concentration readings between the NanoDrop Spectrophotometer and Qubit Fluorometer: **High-quality DNA should show relative agreement in concentration measurements.**
- If you observe a large difference in the concentration readings between the NanoDrop and Qubit systems, for example a difference of **greater than or equal to 50%** of the Qubit dsDNA assay reading, first check for RNA contamination using the Qubit RNA broad range assay.
- If there is no RNA contamination, then we recommend doing at least **one to three rounds of AMPure PB bead purification** until the concentrations are less than 50% different.
- If the agreement does not improve after three rounds of purification, try using either a commercial kit, isopropanol precipitation, or a new DNA extraction method to obtain a cleaner DNA sample.
- If there is RNA in the sample, then treat with **RNase A** followed by a round of AMPure PB bead purification.



Recommended cleanup process for isolated gDNA using 0.45X AMPure PB beads.\* (1X AMPure PB beads may also be used.)



# Methods for evaluation of genomic DNA quality

# Methods for evaluation of genomic DNA quality

Use recommended tools for evaluation of gDNA quality to generate optimal SMRT sequencing data quality

## DNA Sizing QC



Use a Femto Pulse system or pulsed-field gel electrophoresis (PFGE) system for accurate DNA sizing QC of gDNA samples

- ✓ **High-quality, high-molecular weight gDNA** → Longer read lengths / higher data yields
- ✗ **Low-quality, degraded/damaged gDNA** → Shorter read lengths / lower data yields / lower library synthesis yields

## DNA Purity QC



Use a NanoDrop instrument or other spectrophotometer device to determine DNA purity

- ✓ **High-quality, pure gDNA** → Longer read lengths / higher data yields
- ✗ **Low-quality, contaminated DNA** → Shorter read lengths / lower data yields / lower library synthesis yields

## DNA Quantification QC



Use a Qubit fluorometric assay for accurate dsDNA quantitation

- ✓ **Accurate dsDNA quantitation** → Optimal library construction yields / higher data yields
- ✗ **Inaccurate dsDNA quantitation** → Lower library construction yields / lower data yields

# Genomic DNA sizing characterization

## Femto Pulse System [ [Agilent Technologies](#) ]



**Highly Recommended**

Resolves up to ~165 kb  
Requires <1 ng of sample  
<1.5-hour analysis time

## PippinPulse System [ [Sage Science](#) ]




Resolves up to ~80 kb  
Requires ≥50 ng of sample  
~16-hour analysis time


## CHEF Mapper XA System [ [Bio-Rad](#) ]



Resolves up to ~10 Mb  
Requires ≥100 ng of sample  
~16-h analysis time



Lane 1: 8-48 kb Ladder (Bio-Rad)  
Lane 2: 5 kb Ladder (Bio-Rad)  
Lane 3: HMW gDNA  
Lane 4: Degraded gDNA



Lane 1: High MW gDNA  
Lane 2: Degraded gDNA  
Lane 3: 165 kb Ladder

Evaluation of gDNA quality using A) Bio-Rad CHEF Mapper System and B) Femto Pulse System. Lanes A3 and B1 are examples of high quality, high-molecular weight genomic DNA. Lanes A4 and B2 are examples of degraded genomic DNA.

Refer to the appropriate [PacBio Procedure & checklist](#) protocol [documentation](#) for recommended minimum input gDNA fragment size distributions for specific whole genome sequencing applications



# DNA purity determination

- DNA purity can be determined by using a **NanoDrop** system [[Thermo Fisher Scientific](#)] or other spectrophotometer tool
- For ultrapure DNA, A260/280 ratio is typically between ~1.8 - 2.0 and A260/230 ratio is  $\geq 2.0$
- High UV absorbance values are *not* always a guarantee of optimal sequencing performance because not all inhibitors absorb at the wavelengths of 230, 260, and 280 nm.
- Conversely, low UV absorbance values are *not* always a guarantee that non-optimal sequencing performance will be obtained for a sample
  - PacBio has found that gDNA samples with **A260/280 ratios  $\geq 1.7$**  and **A260/230 ratios  $\geq 1.5$**  can still generate excellent HiFi sequencing performance.\*
  - For WGS samples with absorbance ratios outside this range, we recommended performing a **1X AMPure PB bead cleanup** to remove potential contaminants.

## A260/A280 Ratio

- A **low** A260/A280 ratio may be the result of:
  - Protein
  - Phenol
  - Other contaminants that absorb strongly at or near 280 nm
  - Sometimes it may be caused by a very low concentration of nucleic acid.
- **High** 260/280 ratios are not indicative of an issue

If A260/280 and A260/230 readings are out of the recommended ranges, perform one or more rounds of **AMPure PB bead purification** followed by re-assessment of the quantity and purity of the input DNA sample.\*\*\*

## A260/A230 Ratio

- A **low** A260/A230 ratio may be the result of:
  - Protein\*\*
  - Carbohydrate carryover (often a problem with plants)
  - Residual phenol from nucleic acid extraction
  - Residual guanidine (often used in column-based kits)
  - Glycogen used for precipitation
- A **high** A260/A230 ratio may be the result of:
  - Making a blank measurement on a dirty pedestal of a Nanodrop instrument
  - Using an inappropriate solution for the blank measurement



\* See PacBio Technical Note: Sample preparation for PacBio HiFi sequencing from human whole blood ([102-326-500](#))

\*\* See NEB Technical Note: A Practical Guide to Analyzing Nucleic Acid Concentration and Purity with Microvolume Spectrophotometers ([2019](#))

\*\*\* See PacBio Technical Note: Preparing DNA for PacBio HiFi sequencing - Extraction and quality control ([102-193-651](#))

# DNA quantification

- **Accurate** measurement of DNA concentration is critical for PacBio template preparation procedures.
  - Specifically, it is critical to determine the concentration of the **double-stranded DNA**, since only double-stranded DNA will be converted into sequencing templates.
- PacBio highly recommends using a **Qubit fluorometer system** [[Thermo Fisher Scientific](#)] for DNA quantitation
  - Can use a **Qubit dsDNA broad range (BR) assay kit** for initial genomic DNA QC evaluation
  - Use a **Qubit dsDNA high sensitivity (HS) assay kit** for routine DNA quantitation during SMRTbell library construction and sequencing preparation using SMRT Link Sample Setup.
- When assessing gDNA QC, PacBio recommends using **both** fluorometric and spectrophotometric methods – for example, using both the Qubit and NanoDrop instruments
  - If the sample is pure gDNA, free of any RNA contaminants and other small molecules, the two methods should converge to similar DNA concentration measurement values

If the measured NanoDrop concentration is significantly different (>50%) from the Qubit measurement, we recommend doing an AMPure PB bead purification step, followed by a re-measurement with both methods. Typically, a single AMPure PB Bead purification step resolves the discrepancy.

- If the concentration measurement discrepancy after one or more rounds of AMPure PB bead purification is not reduced, we recommend trying another cleanup approach before a re-measurement with both methods.\*





# Storage and shipping of genomic DNA and SMRTbell libraries

# Guidelines for storage and shipping of genomic DNA and SMRTbell libraries

## High-molecular weight genomic DNA storage

- Very clean HMW gDNA can be stored at 4°C for weeks with no degradation. It can also be stored a few days at 25°C.
- HMW gDNA can also be frozen at -20/-80°C for extended storage. **Avoid freeze/thaw cycles.**

## SMRTbell library storage

- If planning to sequence within ~1 week of library generation, storing the SMRTbell library at 4°C in Elution buffer (EB) is recommended
- For storage longer than 1 week, aliquot the SMRTbell library and store at -20°C
- **Reduce or eliminate freeze/thaw cycles** of your SMRTbell library to prevent damage

## Shipping

- **Heat exposure to DNA should be minimized** or eliminated during transport (incubation at 37° C for 1 hour has been shown to cause DNA damage that may result in impaired sequencing performance)
- Lyophilized DNA may be used as long as heat is not applied during the process
- PacBio generally recommends shipping genomic DNA and SMRTbell libraries in a frozen state on dry ice\* with overnight shipping priority
  - Place the primary sample tube(s) inside a secondary form of containment like a 50 mL conical tube and surround it with bubble wrap to help ensure that the primary sample tube does not become damaged during transport

\* **Note:** Genomic DNA extracted with [Circulomics](#) Nanobind Kits can also be shipped in liquid form on wet ice.





# DNA sample extraction literature resources

# DNA sample extraction documentation & other literature

## Technical notes

- Technical note: Preparing samples for PacBio whole genome sequencing for de novo assembly – Collection and storage ([TN100-040518](#))
- Technical Note: Preparing DNA for PacBio HiFi sequencing – Extraction and quality control ([102-193-651](#))
- Technical note: Sample preparation for PacBio HiFi sequencing from human whole blood ([102-326-500](#))

## Nanobind kit handbooks

- Nanobind CBB Big DNA Kit Handbook ([HBK-CBB-001](#))
- Nanobind Plant Nuclei Big DNA Kit Handbook ([HBK-PLT-001](#))
- Nanobind Tissue Big DNA Kit Handbook ([HBK-TIS-001](#))



## Nanobind protocols

- HMW DNA Extraction – Animal Tissue Protocol (Standard Dounce Homogenizer) ([EXT-DHH-001](#))
- HMW DNA Extraction – Animal Tissue Protocol (Standard TissueRuptor) ([EXT-TRH-001](#))
- HMW DNA Extraction – Cultured Cells Protocol ([EXT-CLH-001](#))
- HMW DNA Extraction – Gram Negative Bacteria Protocol ([EXT-GNH-001](#))
- HMW DNA Extraction – Gram Positive Bacteria Protocol ([EXT-GPH-001](#))
- HMW DNA Extraction – Mammalian Whole Blood Protocol (200  $\mu$ L) ([EXT-BLH-001](#))
- HMW DNA Extraction – Nucleated Blood Protocol ([EXT-NBH-001](#))
- HMW DNA Extraction – Plant Nuclei Protocol ([EXT-PLH-001](#))
- Nuclei Isolation – Plant Tissue Protocol (TissueRuptor) ([NUC-TRP-001](#))
- Nuclei Isolation – Plant Tissue Protocol (LN2) ([NUC-LNP-001](#))

# DNA sample extraction documentation & other literature (cont.)

## Nanobind application notes

- cryoPREP Tissue Homogenization Application Note [ [Link](#) ]
- Human Breast Application Note [ [Link](#) ]
- Mammalian Brain Application [ [Link](#) ]
- Mammalian Liver Application [ [Link](#) ]
- Mammalian Spleen Application Note [ [Link](#) ]
- Fish Testis Application Note [ [Link](#) ]
- Snail Application Note [ [Link](#) ]
- Crab Application Note [ [Link](#) ]
- Fish Skeletal Muscle Application Note [ [Link](#) ]
- Aplysia Application Note [ [Link](#) ]

### Have challenging or interesting samples to process?

Please contact [Circulomics](#) for recommendations on how to process your challenging or interesting samples using Nanobind Kits.

## Nanobind automation scripts [contact [Circulomics](#)]

- Nanobind KF CBB Big DNA Kit Automated Reference Guide
- Automated HMW DNA Extraction – Cultured Cells
- Automated HMW DNA Extraction – 200 µL Mammalian Blood
- Automated HMW DNA Extraction – 1 mL Mammalian Blood
- Automated HMW DNA Extraction – Gram-Negative Bacteria
- Automated HMW DNA Extraction – Gram-Positive Bacteria
- Automated HMW DNA Extraction – Animal Tissue



Circulomics Technical Support: [support@circulomics.com](mailto:support@circulomics.com)

# DNA sample preparation online resource

Literature resource for sample collection and DNA extraction protocol references

The listing below is a collection of publications by the scientific community describing extraction protocols for high-molecular weight DNA followed by PacBio sequencing. When possible, the links point directly to the methods section (or supplementary information).

[Animals](#)   [Plants](#)   [Fungi](#)   [Protists](#)

If you have protocols you would like to share, or have questions about DNA extraction for PacBio sequencing, please contact [ExtractDNA@pacb.com](mailto:ExtractDNA@pacb.com).

## Animals

### Invertebrates

- [Panova2016](#) – DNA extraction protocols for whole-genome sequencing in marine organisms
- [microinvertebrates](#)
  - [Laumer2020](#) – protocols for diverse meiofauna species, including *C. elegans* & *Caecella truncata* (SMRT Leiden conference presentation)
- [arthropods](#)
  - [arachnids](#)
    - [Guerrero2019](#) – The Pacific Biosciences de novo assembled genome dataset from a parthenogenetic New Zealand wild population of the longhorned tick, *Haemaphysalis longicornis* Neumann, 1901



### Crustaceans

- [Liu2019](#) – Crab muscle extraction & sequencing (Circulomics application note)
- [Baldwin-Brown2018](#) – A New Standard for Crustacean Genomes: The Highly Contiguous, Annotated Genome Assembly of the Clam Shrimp *Enlimnadia texana* Reveals HOX Gene Order and Identifies the Sex Chromosome
- [Gonçalves Athanasio2016](#) – Optimisation of DNA extraction from the crustacean *Daphnia*

## Plants

### Methods

#### Sample collection

A female yellowbelly pufferfish (Fig. 2), reared in the fish breeding centre of Fujian Normal University in Fuzhou City of Fujian Province was used for genome sequencing and assembly. Fresh white muscle, eye, skin, gonad, gut, liver, kidney, blood, gall bladder and air bladder tissues were collected and quickly frozen in liquid nitrogen for one hour. White muscle tissues were used for DNA sequencing for genome assembly, while all tissues were used for transcriptome sequencing.

Fig. 2



A picture of the yellowbelly pufferfish used in the genome sequencing and assembly.

#### DNA and RNA sequencing

Genomic DNA from white muscle tissue was extracted using the

at genomic DNA for long-read sequencing of single molecules

ce Genome With Single-Molecule Sequencing Uncovers a Recent Burst of Genes

sembly reveals the structure of the *Arabidopsis thaliana* Nd-1 genome and its gene set

sive, high-throughput plant genomic DNA extraction method suitable for

or consist

ing and

ect run

olly pu

ace for the Fast-Growing Microalga *Picochlorum celeri*

*Undaria pinnatifida*: Chromosome-Level Assembly Using PacBio and Hi-C

ne suggests convergent functions of homeobox genes in algae and land plants (algae)

whole-genome sequencing in marine organisms (algae)

[www.ExtractDNAforPacBio.com](http://www.ExtractDNAforPacBio.com)

## Fungi

- [Faure2019](#) – Long-Read Genome Sequence of the Sugar Beet Rhizosphere Mycoparasite *Pythium oligandrum*
- [Nagappan2018](#) – Improved nucleic acid extraction protocols for *Candida boninense*, *G. miniatacinctum* and *G. torulosa*
- [Schwessinger2017](#) – Extraction of high molecular weight DNA from fungal rust spores for long read sequencing
- [Solomon2016](#) – Robust and effective methodologies for cryopreservation and DNA extraction from anaerobic gut fungi
- [Sonnenberg2016](#) – A detailed analysis of the recombination landscape of the button mushroom *Agaricus bisporus* var. *bisporus*



# PacBio

[www.pacb.com](http://www.pacb.com)

Research use only. Not for use in diagnostic procedures. © 2022, PacBio. All rights reserved. Pacific Biosciences, the PacBio logo, PacBio, Circulomics, Omnione, SMRT, SMRTbell, Iso-Seq, Sequel, Nanobind, and SBB are trademarks of Pacific Biosciences of California, Inc. (PacBio). All other trademarks are the sole property of their respective owners.