

3M Science.
Applied to Life.™

Advancing food safety through smart innovation.

3M Food Safety
Full Product Catalog

June 2021

Table of Contents


3M™ Petrifilm™ Plates — Indicator Organism Testing	2–4
3M™ Petrifilm™ Standard Plates	
3M™ Petrifilm™ Rapid Plates	
3M™ Petrifilm™ Spreaders	
3M™ Petrifilm™ Plate Reader Advanced	
3M™ Molecular Detection System for Pathogen Testing — Pathogen Detection	5–6
3M™ Molecular Detection System Assays	
3M™ Molecular Detection System Accessories	
3M™ Molecular Detection System Enrichment Media	
3M™ Petrifilm™ <i>Salmonella</i> Express System — Pathogen Detection	7
3M™ Petrifilm™ <i>Salmonella</i> Enrichment Base and <i>Salmonella</i> Enrichment Supplement	
3M™ Clean-Trace™ Hygiene Monitoring and Management System — Hygiene Monitoring	8–9
3M™ Clean-Trace™ Hygiene Monitoring and Management System LM1	
3M™ Clean-Trace™ Surface ATP Tests	
3M™ Clean-Trace™ Water ATP Tests	
3M™ Clean-Trace™ Protein Tests	
3M™ Allergen Protein Tests — Allergen Detection	10–11
3M™ Allergen Protein Rapid Kits	
3M™ Allergen Protein Enzyme-Linked Immunosorbent Assay (ELISA) Kits	
3M™ Allergen Protein Testing Accessories	
3M™ Sample Collection, Preparation & Processing	12–17
Collect	
Scrub Samplers	
Sponges	
Swabs	
Plain Sample Bags	
Printed Sample Bags	
Carcass Sampling	
Prepare	
Flip-Top Dilution Bottles	
Mini-Flip Top Dilution Vials	
Blender Bags	
Dilution and Enrichment Pouches	
Dehydrated Media	
Process	
General Lab Equipment	
3M™ Microbial Luminescence System	18–19
3M™ Microbial Luminescence Systems	
3M™ Microbial Luminescence System Accessories	
3M™ Time & Temperature Indicators	20
3M™ Freeze Watch™ Indicators	
3M™ MonitorMark™ Time Temperature Indicators	
3M™ Food Safety Hardware Service	21

Indicator Organism Testing

3M™ Petrifilm™ Plates

Fresh science for the modern lab.

Recognized with over 100 global recognitions and validations, 3M™ Petrifilm™ Plates are ready-to-use plates that reduce or eliminate the time-consuming agar prep process. With proven testing methods for consistent, reliable results, these compact plates can increase productivity, use less storage/incubator space and generate less waste.



3M™ Petrifilm™ Standard Plates



3M™ Petrifilm™ Aerobic Count Plate



3M™ Petrifilm™ Coliform Count Plate



3M™ Petrifilm™ *E. coli*/Coliform Count Plate



3M™ Petrifilm™ *Enterobacteriaceae* Count Plate



3M™ Petrifilm™ Environmental *Listeria* Plate



3M™ Petrifilm™ High-Sensitivity Coliform Count Plate



3M™ Petrifilm™ Lactic Acid Bacteria Count Plate



3M™ Petrifilm™ Select *E. coli* Count Plate



3M™ Petrifilm™ Staph Express Count Plate and Disk



3M™ Petrifilm™ Yeast & Mold Count Plate



3M™ Petrifilm™ Aqua Heterotrophic Count Plate

3M™ Petrifilm™ Rapid Plates



3M™ Petrifilm™ Rapid
Aerobic Count Plate



3M™ Petrifilm™ Rapid
Coliform Count Plate



3M™ Petrifilm™ Rapid
E. coli/Coliform Count Plate



3M™ Petrifilm™ Rapid
Yeast & Mold Count Plate

3M™ Petrifilm™ Spreaders



3M™ Petrifilm™
Environmental *Listeria*
Plate Spreader



3M™ Petrifilm™
High-Sensitivity
Plate Spreader



3M™ Petrifilm™
Flat Spreader

The 3M™ Petrifilm™ Spreaders are easy-to-use tools designed to quickly and evenly distribute inoculum within the 3M™ Petrifilm™ Plates. The spreaders are reusable and do not contact plating surface or sample.

3M™ Petrifilm™ Plate Reader Advanced

Automatically reads and records results of 3M™ Petrifilm™ Plates in four to six seconds. Automatic data storage saves time and mitigates transcription errors.



Quality Indicator Testing — 3M™ Petrifilm™ Plates

Catalog Number	Description	Quantity/Billing Unit
3M™ Petrifilm™ Standard Plates		
6400	3M™ Petrifilm™ Aerobic Count Plates 6400	100/box
6406	3M™ Petrifilm™ Aerobic Count Plates 6406	1,000/case
6450	3M™ Petrifilm™ Aqua Heterotrophic Count Plates 6450	100/box
6452	3M™ Petrifilm™ Aqua Heterotrophic Count Plates 6452	1,000/case
6410	3M™ Petrifilm™ Coliform Count Plates 6410	50/box
6416	3M™ Petrifilm™ Coliform Count Plates 6416	1,000/case
6404	3M™ Petrifilm™ <i>E. coli</i> /Coliform Count Plates 6404	50/box
6414	3M™ Petrifilm™ <i>E. coli</i> /Coliform Count Plates 6414	500/case
6420	3M™ Petrifilm™ <i>Enterobacteriaceae</i> Count Plates 6420	50/box
6421	3M™ Petrifilm™ <i>Enterobacteriaceae</i> Count Plates 6421	1,000/case
6447	3M™ Petrifilm™ Environmental <i>Listeria</i> Plates 6447	50/box
6448	3M™ Petrifilm™ Environmental <i>Listeria</i> Plates 6448	200/case
6405	3M™ Petrifilm™ High-Sensitivity Coliform Count Plates 6405	50/box
6415	3M™ Petrifilm™ High-Sensitivity Coliform Count Plates 6415	500/case
6461	3M™ Petrifilm™ Lactic Acid Bacteria Count Plates 6461	50/box
6462	3M™ Petrifilm™ Lactic Acid Bacteria Count Plates 6462	500/case
6434	3M™ Petrifilm™ Select <i>E. coli</i> Count Plates 6434	50/box
6435	3M™ Petrifilm™ Select <i>E. coli</i> Count Plates 6435	500/case
6490	3M™ Petrifilm™ Staph Express Count Plates 6490	50/box
6491	3M™ Petrifilm™ Staph Express Count Plate 6491	500/case
6492	3M™ Petrifilm™ Staph Express Disks 6492	20/box
6493	3M™ Petrifilm™ Staph Express Disks 6493	100/case
6407	3M™ Petrifilm™ Yeast & Mold Count Plates 6407	100/box
6417	3M™ Petrifilm™ Yeast & Mold Count Plates 6417	1,000/case
3M™ Petrifilm™ Rapid Plates		
6478	3M™ Petrifilm™ Rapid Aerobic Count Plates 6478	50/box
6479	3M™ Petrifilm™ Rapid Aerobic Count Plates 6479	500/case
6402	3M™ Petrifilm™ Rapid Coliform Count Plates 6402	50/box
6412	3M™ Petrifilm™ Rapid Coliform Count Plates 6412	500/case
6436	3M™ Petrifilm™ Rapid <i>E. coli</i> /Coliform Count Plate 6436	50/box
6437	3M™ Petrifilm™ Rapid <i>E. coli</i> /Coliform Count Plate 6437	500/case
6475	3M™ Petrifilm™ Rapid Yeast & Mold Count Plates 6475	50/box
6477	3M™ Petrifilm™ Rapid Yeast & Mold Count Plates 6477	500/case
3M™ Petrifilm™ Spreaders		
6498	3M™ Petrifilm™ Environmental <i>Listeria</i> Plate Spreader 6498	2/case
6481	3M™ Petrifilm™ High-Sensitivity Plate Spreader 6481	2/case
6425	3M™ Petrifilm™ Flat Spreader 6425	2/case
3M™ Petrifilm™ Plate Reader Advanced		
6558	3M™ Petrifilm™ Plate Reader Advanced 6558	1 kit/case

Pathogen Detection

3M™ Molecular Detection System for Pathogen Testing

Make critical decisions faster with the 3M Molecular Detection System.

3M™ Molecular Detection System for Pathogen Testing offers a broad range of pathogen testing to detect foodborne pathogens in a wide variety of foods and food environments. Each compact, streamlined component works together to optimize performance with features that save time and increase productivity.



3M™ Molecular Detection System Assays



3M™ Molecular Detection Assay 2 – *Salmonella*



3M™ Molecular Detection Assay 2 – *Listeria monocytogenes*



3M™ Molecular Detection Assay 2 – *Listeria*



3M™ Molecular Detection Assay 2 – *Campylobacter* with 3M™ *Campylobacter* Enrichment Broth CE250



3M™ Molecular Detection Assay 2 – *E. coli* O157 (including H7)



3M™ Molecular Detection Assay 2 – STEC Gene Screen (stx and eae)



3M™ Molecular Detection Assay 2 – STEC Gene Screen (stx)



3M™ Molecular Detection Assay 2 – *Cronobacter*

3M™ Molecular Detection System Accessories




3M™ Molecular Detection Matrix Control




3M™ Molecular Detection Cap/Decap Tool – Lysis



3M™ Molecular Detection Cap/Decap Tool – Reagent



3M™ Molecular Detection Chill Block Insert



3M™ Molecular Detection Heat Block Insert



3M™ Molecular Detection Speed Loader Tray

Pathogen Detection — 3M™ Molecular Detection System for Pathogen Testing

Catalog Number	Description	Quantity/Billing Unit
3M™ Molecular Detection System Assays		
MDA2CAM96	3M™ Molecular Detection Assay 2 – <i>Campylobacter</i> MDA2CAM96	96 tests/kit
MDA2CRO96	3M™ Molecular Detection Assay 2 – <i>Cronobacter</i> MDA2CRO96	96 tests/kit
MDA2ECO96	3M™ Molecular Detection Assay 2 – <i>E. coli</i> O157 (including H7) MDA2ECO96	96 tests/kit
MDA2LIS96	3M™ Molecular Detection Assay 2 – <i>Listeria</i> MDA2LIS96	96 tests/kit
MDA2LMO96	3M™ Molecular Detection Assay 2 – <i>Listeria monocytogenes</i> MDA2LMO96	96 tests/kit
MDA2SAL96	3M™ Molecular Detection Assay 2 – <i>Salmonella</i> MDA2SAL96	96 tests/kit
MDA2STXEAE48	3M™ Molecular Detection Assay 2 – STEC Gene Screen (<i>stx</i> and <i>eae</i>), MDA2STXEAE48	48 tests/kit
MDA2STX96	3M™ Molecular Detection Assay 2 – STEC Gene Screen (<i>stx</i>), MDA2STX96	96 tests/kit
3M™ Molecular Detection System Accessories		
MDSCDL	3M™ Molecular Detection Cap/Decap Tool – Lysis MDSCDL	1/case
MDSCDR	3M™ Molecular Detection Cap/Decap Tool – Reagent MDSCDR	1/case
MDSCBIN	3M™ Molecular Detection Chill Block Insert MDSCBIN	1/case
MDSHBIN	3M™ Molecular Detection Heat Block Insert MDSHBIN	1/case
MDS100	3M™ Molecular Detection Instrument and Accessory Kit MDS100	1/case
MDMC96	3M™ Molecular Detection Matrix Control MDMC96	96 tests/kit
MDSSLT	3M™ Molecular Detection Speed Loader Tray MDSSLT	1/case
3M™ Molecular Detection System Enrichment Media		
CE250	3M™ <i>Campylobacter</i> Enrichment Broth CE250, 250g	1/case

Pathogen Detection

3M™ Petrifilm™ *Salmonella* Express System

Faster, simpler, all-in-one method for *Salmonella* pathogen testing that uses 3M™ Petrifilm™ Plate technology. Plate provides presumptive positive *Salmonella* results in as little as 40 hours for low-background, processed foods and 48 hours for high-background, non-processed foods. An all-in-one test and biochemical confirmation method. Confirmation disk easily added to plate for biochemical confirmation provides results in 4–5 hours.



3M™ Petrifilm™ *Salmonella* Enrichment Base and *Salmonella* Enrichment Supplement



3M™ *Salmonella* Enrichment Base



3M™ *Salmonella* Enrichment Supplement

Pathogen Detection — 3M™ Petrifilm™ *Salmonella* Express System

Catalog Number	Description	Quantity/Billing Unit
3M™ Petrifilm™ <i>Salmonella</i> Express Systems		
6536	3M™ Petrifilm™ <i>Salmonella</i> Express Plate 6536	50/box
6537	3M™ Petrifilm™ <i>Salmonella</i> Express Plate 6537	200/case
6538	3M™ Petrifilm™ <i>Salmonella</i> Express Confirmation Disk 6538	5/box
6539	3M™ Petrifilm™ <i>Salmonella</i> Express Confirmation Disk 6539	25/case
SEB500	3M™ <i>Salmonella</i> Enrichment Base SEB500, 500g	1/case
SESUP001	3M™ <i>Salmonella</i> Enrichment Supplement SESUP001, 1g	1/case

Hygiene Monitoring

3M™ Clean-Trace™ Hygiene Monitoring and Management System

See beyond the surface to evaluate the effectiveness of your cleaning procedures.

The 3M™ Clean-Trace™ Hygiene Monitoring and Management System is a robust detection system that quickly provides stable, consistent results and helps you leverage data. The 3M™ Clean-Trace™ Software turns hygiene monitoring data into actionable information to help you find high-risk areas in your food processing facilities before you begin production. The ergonomic design of the 3M™ Clean-Trace™ Luminometer makes testing simple and can help to minimize training time and costs.



Hygiene Monitoring — 3M™ Clean-Trace™ Hygiene Monitoring and Management System

Catalog Number	Description	Quantity/Billing Unit
3M™ Clean-Trace™ Hygiene Monitoring and Management System		
LM1	3M™ Clean-Trace™ Hygiene Monitoring and Management System LM1	1/case

3M™ Clean-Trace™ Surface ATP Tests

Help detect the presence of adenosine triphosphate (ATP) to verify cleaning effectiveness. 3M™ Clean-Trace™ Surface ATP Tests are easy to use and implement. The sleek, portable swabs are easy to activate with one hand.




3M™ Clean-Trace™ Surface ATP

3M™ Clean-Trace™ Water ATP Tests

Help detect the presence of adenosine triphosphate (ATP) in rinse water to verify cleaning effectiveness to verify cleaning effectiveness. 3M™ Clean-Trace™ Water ATP Tests are easy to use and implement. The sleek, portable swabs are easy to activate with one hand.



3M™ Clean-Trace™ Water Plus – Total ATP




3M™ Clean-Trace™ Water – Free ATP

3M™ Clean-Trace™ Protein Tests

Help detect the presence of protein on surfaces and in solution. Detect as little as 3µg of protein. Test for the presence of protein to verify effectiveness of cleaning; visual reading of a color change indicates the presence of detectable levels of protein.



3M™ Clean-Trace™ Surface Protein Plus



3M™ Clean-Trace™ Surface Protein (Allergen)

Hygiene Monitoring — 3M™ Clean-Trace™ Hygiene Monitoring and Management System (*cont.*)

Catalog Number	Description	Quantity/Billing Unit
3M™ Clean-Trace™ Surface ATP Tests		
UXL100	3M™ Clean-Trace™ Surface ATP Test Swab UXL100	100/case
ATP10	3M™ Clean-Trace™ Surface Positive Control ATP10	10/case
3M™ Clean-Trace™ Water ATP Tests		
AQF100	3M™ Clean-Trace™ Water – Free ATP AQF100	100/case
AQT200	3M™ Clean-Trace™ Water Plus – Total ATP AQT200	100/case
HWATP10	3M™ Clean-Trace™ Water Positive Control HWATP10	10/case
LWATP10	3M™ Clean-Trace™ Water Positive Control LWATP10	10/case
3M™ Clean-Trace™ Protein Tests		
PRO100	3M™ Clean-Trace™ Surface Protein Plus Test Swab PRO100	100/case
ALLTEC60	3M™ Clean-Trace™ Surface Protein (Allergen) Test Swab ALLTEC60	60/case

Allergen Detection

3M™ Allergen Protein Tests

Complex testing made simple.

3M™ Allergen Testing products are designed to detect processed and unprocessed target allergen proteins and allow for a simplified workflow. Tests can accommodate environmental swab, CIP rinse water and food samples.



*Lateral Flow Devices
for Specific Allergen
Detection*

3M™ Allergen Protein Rapid Tests

Rapid, qualitative immunochromatographic lateral flow test cartridges designed for the detection of proteins from allergenic foods.



3M™ Allergen Protein Enzyme-Linked Immunosorbent Assay (ELISA) Tests

Confirmed, quantifiable results with solid phase sandwich ELISA tests for monitoring specific food allergens across a wide assortment of foods.



Allergen Detection — 3M™ Allergen Protein Tests

Catalog Number	Description	Quantity/Billing Unit
3M™ Allergen Protein Rapid Kits		
L25ALM	3M™ Almond Protein Rapid Kit L25ALM	25 tests/kit
L25CHW	3M™ Cashew Protein Rapid Kit L25CHW	25 tests/kit
L25COC	3M™ Coconut Protein Rapid Kit L25COC	25 tests/kit
L25EGG	3M™ Egg White Protein Rapid Kit L25EGG	25 tests/kit
L25FSH	3M™ Fish Protein Rapid Kit L25FSH	25 tests/kit
L25GLU	3M™ Gluten Protein Rapid Kit L25GLU	25 tests/kit
L25HZL	3M™ Hazelnut Protein Rapid Kit L25HZL	25 tests/kit
L25MLK	3M™ Milk Protein Rapid Kit L25MLK	25 tests/kit
L25PNT	3M™ Peanut Protein Rapid Kit L25PNT	25 tests/kit
L25PEC	3M™ Pecan Protein Rapid Kit L25PEC	25 tests/kit
L25PST	3M™ Pistachio Protein Rapid Kit L25PST	25 tests/kit
L25SOY	3M™ Soy Protein Rapid Kit L25SOY	25 tests/kit
L25WAL	3M™ Walnut Protein Rapid Kit L25WAL	25 tests/kit
3M™ Allergen Protein ELISA Kits		
E96ALM	3M™ Almond Protein ELISA Kit E96ALM	96 wells/kit
E96CAS	3M™ Bovine Casein Protein ELISA Kit E96CAS	96 wells/kit
E96MLK	3M™ Bovine Total Milk Protein ELISA Kit E96MLK	96 wells/kit
E96BZL	3M™ Brazil Nut Protein ELISA Kit E96BZL	96 wells/kit
E96CHW	3M™ Cashew Protein ELISA Kit E96CHW	96 wells/kit
E96COC	3M™ Coconut Protein ELISA Kit E96COC	96 wells/kit
E96CRU	3M™ Crustacean Protein ELISA Kit E96CRU	96 wells/kit
E96EGG	3M™ Egg White Protein ELISA Kit E96EGG	96 wells/kit
E96FSH	3M™ Fish Protein ELISA Kit E96FSH	96 wells/kit
E96HZL	3M™ Hazelnut Protein ELISA Kit E96HZL	96 wells/kit
E96MAC	3M™ Macadamia Protein ELISA Kit E96MAC	96 wells/kit
E96MOL	3M™ Mollusk Protein ELISA Kit E96MOL	96 wells/kit
E96MUS	3M™ Mustard Protein ELISA Kit E96MUS	96 wells/kit
E96PNT	3M™ Peanut Protein ELISA Kit E96PNT	96 wells/kit
E96PEC	3M™ Pecan Protein ELISA Kit E96PEC	96 wells/kit
E96PNE	3M™ Pine Nut Protein ELISA Kit E96PNE	96 wells/kit
E96PST	3M™ Pistachio Protein ELISA Kit E96PST	96 wells/kit
E96SES	3M™ Sesame Protein ELISA Kit E96SES	96 wells/kit
E96SOY	3M™ Soy Protein ELISA Kit E96SOY	96 wells/kit
E96WAL	3M™ Walnut Protein ELISA Kit E96WAL	96 wells/kit
3M™ Allergen Protein Testing Accessories		
ECS100	3M™ Environmental Collection Swab ECS100	100/case
TP100	3M™ Transfer Pipette TP100	100/pack

3M™ Sample Collection, Preparation & Processing

Consistency from the start.

Better testing begins with better samples — and better samples are obtained when using quality sample collection tools. 3M™ Sample Collection provides quality products designed specifically for microbial surface sampling.



Collect. Sponges and swabs designed specifically for microbial surface sampling.



Transport. A wide variety of bags for use during transport to testing lab.



Prepare. Standardize testing results with a wide variety of prepared media, in convenient flip-top vials, pouches or bottles.

Sample Collection Products



3M™ Environmental Scrub Sampler with 10mL Wide Spectrum Neutralizer




3M™ Hydrated-Sponges and 3M™ Sponge-Sticks



3M™ Quick Swab



3M™ Enviro Swab



3M™ Swab-Samplers



3M™ Sample Bags



3M™ Carcass Sample Kit

Sample Preparation & Processing Products



3M™ Flip-Top Dilution Bottles



3M™ Mini Flip-Top Dilution Vials



3M™ Blender Bags



3M™ Dilution Pouches and 3M™ Enrichment Pouches



3M™ Dehydrated Media

Sample Collection, Preparation & Processing — Collect

Catalog Number	Description	Quantity/Billing Unit
Scrub Samplers		
ESS10WSN	3M™ Environmental Scrub Sampler Stick with 10mL Wide Spectrum Neutralizer	100/Case
HES10WSN2G	3M™ Environmental Scrub Sampler with 10mL Wide Spectrum Neutralizer and Gloves	100/Case
Sponges		
BP133ES	3M™ Dry-Sponge BP133ES	100/case
BP237SPG	3M™ Dry-Sponge with Glove BP237SPG	100/case
HS10BPW	3M™ Hydrated-Sponge with 10mL Buffered Peptone Water Broth HS10BPW	100/case
HS10BPW2G	3M™ Hydrated-Sponge with 10mL Buffered Peptone Water Broth and Gloves HS10BPW2G	100/case
HS10DE2G	3M™ Hydrated-Sponge with 10mL D/E Neutralizing Buffer and Gloves HS10DE2G	100/case
HS10LET	3M™ Hydrated-Sponge with 10mL Lethen Broth HS10LET	100/case
HS10LET2G	3M™ Hydrated-Sponge with 10mL Lethen Broth and Gloves HS10LET2G	100/case
HS10NB	3M™ Hydrated-Sponge with 10mL Neutralizing Buffer HS10NB	100/case
HS10NB2G	3M™ Hydrated-Sponge with 10mL Neutralizing Buffer and Gloves HS10NB2G	100/case
HS2410NB2G	3M™ Hydrated-Sponge with 10mL Neutralizing Buffer and Gloves HS2410NB2G	100/case
SSL100	3M™ Sponge-Stick SSL100	100/case
SSL10BPW	3M™ Sponge-Stick with 10mL Buffered Peptone Water Broth SSL10BPW	100/case
SSL10DE	3M™ Sponge-Stick with 10mL D/E Neutralizing Broth SSL10DE	100/case
SSL10LET	3M™ Sponge-Stick with 10mL Lethen Broth SSL10LET	100/case
SSL10LET2G	3M™ Sponge-Stick with 10mL Lethen Broth and Gloves SSL10LET2G	100/case
SSL10NB	3M™ Sponge-Stick with 10mL Neutralizing Buffer SSL10NB	100/case
SSL10NB2G	3M™ Sponge-Stick with 10mL Neutralizing Buffer and Gloves SSL10NB2G	100/case
XSL008	3M™ Extendable Pole XSL008, 8 ft.	1/case
Swabs		
6433	3M™ Quick Swab 6433	250/box
6432	3M™ Quick Swab 6432	50/box
ENVSWB100	3M™ Enviro Swab ENVSWB100	100/case
RS9601LET	3M™ Swab-Sampler with 1mL Lethen Broth RS9601LET	100/case
RS9604LET	3M™ Swab-Sampler with 4mL Lethen Broth RS9604LET	100/case
RS9605LET	3M™ Swab-Sampler with 5mL Lethen Broth RS9605LET	100/case
RS96010LET	3M™ Swab-Sampler with 10mL Lethen Broth RS96010LET	100/case
RS96010BPW	3M™ Swab-Sampler with 10mL Buffered Peptone Water Broth RS96010BPW	100/case
RS9604NB	3M™ Swab-Sampler with 4mL Neutralizing Buffer RS9604NB	100/case
RS96010NB	3M™ Swab-Sampler with 10mL Neutralizing Buffer RS96010NB	100/case
RS96010DE	3M™ Swab-Sampler with 10mL D/E Neutralizing Broth RS96010DE	100/case

Sample Collection, Preparation & Processing — Collect (cont.)

Catalog Number	Description	Quantity/Billing Unit
Plain Sample Bags		
BP111S	3M™ Plain Sample Bag BP111S, 3 in. x 7 in., 4 oz., 2.5 mil	1,000/case
BP112S	3M™ Plain Sample Bag BP112S, 4.5 in. x 9 in., 15 oz., 2.5 mil	1,000/case
BP113S	3M™ Plain Sample Bag BP113S, 5.5 in. x 9 in., 22 oz., 3 mil	1,000/case
BP118S	3M™ Plain Sample Bag BP118S, 4.5 in. x 12 in., 27 oz., 2.5 mil	1,000/case
BP119S	3M™ Plain Sample Bag BP119S, 4.5 in. x 15 in., 36 oz., 2.5 mil	1,000/case
BP121S	3M™ Plain Sample Bag BP121S, 5.5 in. x 15 in., 50 oz., 3 mil	1,000/case
BP134S	3M™ Plain Sample Bag BP134S, 4.5 in. x 7.5 in., 12 oz., 2.5 mil	1,000/case
BP1015LP	3M™ Plain Sample Bag BP1015LP, 10 in. x 15 in., 135 oz., 4 mil	1,000/case
BP41212	3M™ Plain Sample Bag BP41212, 12 in. x 12 in., 102 oz., 4 mil	250/case
BP41218	3M™ Plain Sample Bag BP41218, 12 in. x 18 in., 228 oz., 4 mil	250/case
BRB3500	3M™ Bird Rinse Plain Sample Bag BRB3500, 15 in. x 20 in., 405 oz., 3 mil	500/case
Printed Sample Bags		
BP114S	3M™ Printable Sample Bag BP114S, 3 in. x 7 in., 4 oz., 2.5 mil	1,000/case
BP115S	3M™ Printable Sample Bag BP115S, 4.5 in. x 9 in., 15 oz., 2.5 mil	1,000/case
BP123S	3M™ Printable Sample Bag BP123S, 3 in. x 5 in., 2 oz., 2.5 mil	1,000/case
BP127S	3M™ Printable Sample Bag BP127S, 5.5 in. x 9 in., 22 oz., 3 mil	1,000/case
BP701	3M™ Printable Sample Bag BP701, 7 in. x 12 in., 55 oz., 3 mil	1,000/case
BP1015	3M™ Printable Sample Bag BP1015, 10 in. x 15 in., 135 oz., 4 mil	1,000/case
Carcass Sampling		
USDA050	3M™ USDA Poultry Template USDA050	100/case
USDA100	3M™ USDA Cattle Template USDA100	100/case
USDA25BFD	3M™ USDA Cattle Screw-top Bottle with 25mL Butterfield's Buffer USDA25BFD	100/case
USDA25BPW	3M™ USDA Cattle Screw-top Bottle with 25mL Buffered Peptone Water Broth USDA25BPW	100/case
USDAKIT025	3M™ USDA Cattle/Swine Kit with 25mL Butterfield's Buffer USDAKIT025	15/case

Sample Collection, Preparation & Processing — Prepare

Catalog Number	Description	Quantity/Billing Unit
Flip-Top Dilution Bottles		
FTBPWI9060	3M™ Flip-Top Dilution Bottle with Buffered Peptone Water (ISO) FTBPWI9060	60/case
FTBFD9060	3M™ Flip-Top Dilution Bottle with Butterfield's Buffer FTBFD9060	60/case
FTBFD9960	3M™ Flip-Top Dilution Bottle with Butterfield's Buffer FTBFD9960	60/case
FTPHB9060	3M™ Flip-Top Dilution Bottle with Phosphate Buffer FTPHB9060	60/case
FTPHB9960	3M™ Flip-Top Dilution Bottle with Phosphate Buffer FTPHB9960	60/case
FTBPW9060	3M™ Flip-Top Dilution Bottle with Buffered Peptone Water FTBPW9060	60/case
FTPW9060	3M™ Flip-Top Dilution Bottle with Peptone Water FTPW9060	60/case
FTPW9960	3M™ Flip-Top Dilution Bottle with Peptone Water FTPW9960	60/case
Mini Flip-Top Dilution Vials		
BPPFV9BPW	3M™ Mini Flip-Top Vial with 9mL Buffered Peptone Water Broth BPPFV9BPW	100/case
BPPFV10BPW	3M™ Mini Flip-Top Vial with 10mL Buffered Peptone Water Broth BPPFV10BPW	100/case
BPPFV9BFD	3M™ Mini Flip-Top Vial with 9mL Butterfield's Buffer BPPFV9BFD	100/case
BPPFV9LB	3M™ Mini Flip-Top Vial with 9mL Lethen Broth BPPFV9LB	100/case
BPPFV10LB	3M™ Mini Flip-Top Vial with 10mL Lethen Broth BPPFV10LB	100/case
BPPFV9MRD	3M™ Mini Flip-Top Vial with 9mL Maximum Recovery Diluent BPPFV9MRD	100/case
BPPFV9MRS	3M™ Mini Flip-Top Vial with 9mL MRS Broth BPPFV9MRS	100/case
BPPFV10DE	3M™ Mini Flip-Top Vial with 10mL D/E Neutralizing Broth BPPFV10DE	100/case
BPPFV10NB	3M™ Mini Flip-Top Vial with 10mL Neutralizing Buffer BPPFV10NB	100/case
Blender Bags		
BP1520F	3M™ Blender Bag with filter (no wire tie) BP1520F, 15 in. x 20 in., 405 oz., 3 mil	200/case
BP912S	3M™ Blender Bag (no wire tie) BP912S, 7 in. x 12 in., 55 oz., 3 mil	1,000/case
BP913S	3M™ Blender Bag (no wire tie) BP913S, 7 in. x 12 in., 55 oz., 4 mil	1,000/case
BP914S	3M™ Blender Bag with Filter (no wire tie) BP914S, 7 in. x 12 in., 55 oz., 3 mil	400/case
6469	3M™ Homogenizer Bag Filtered 6469, 6.75 in. x 12 in., 55 oz., 3 mil	200/case
Dilution and Enrichment Pouches		
QDBFD225F	3M™ Dilution Pouch with 225mL Butterfield's Buffer and Filter QDBFD225F	32/case
QDBFD225	3M™ Dilution Pouch with 225mL Butterfield's Buffer QDBFD225	32/case
QDBFD450	3M™ Dilution Pouch with 450mL Butterfield's Buffer QDBFD450	16/case
QDBFD99	3M™ Dilution Pouch with 99mL Butterfield's Buffer QDBFD99	50/case
QEBPW225	3M™ Enrichment Pouch with 225mL Buffered Peptone Water Broth QEBPW225	32/case
QEBPW225F	3M™ Enrichment Pouch with 225mL Buffered Peptone Water Broth with Filter QEBPW225F	32/case
QEBPWI225F	3M™ Enrichment Pouch with 225mL Buffered Peptone Water Broth ISO with Filter QEBPWI225F	32/case
QEBPWS225	3M™ Enrichment Pouch with 225mL Buffered Peptone Water Broth ISO QEBPWS225	32/case
QELAC225	3M™ Enrichment Pouch with 225mL Lactose Broth QELAC225	32/case
QETSB225	3M™ Enrichment Pouch with 225mL Tryptic Soy Broth QETSB225	32/case

Sample Collection, Preparation & Processing — Prepare (cont.)

Catalog Number	Description	Quantity/Billing Unit
Dehydrated Media		
BP0115005	3M™ Buffered Peptone Water Broth BP0115005, 2.5kg	1/case
BP0115500	3M™ Buffered Peptone Water Broth BP0115500, 500g	1/case
BP0118005	3M™ Demi Fraser Broth Base BP0118005, 2.5kg	4/case
BP0118500	3M™ Demi Fraser Broth Base BP0118500, 500g	1/case
BP0210005	3M™ Fraser Broth Base BP0210005, 2.5kg	4/case
BP0210500	3M™ Fraser Broth Base BP0210500, 500g	1/case
BP0220010	3M™ Fraser Broth Supplement BP0220010, 10mL	10/case
BP0235500	3M™ Lactose Broth BP0235500, 500g	1/case
BP0245500	3M™ Letheen Broth BP0245500, 500g	12/case
BP0255500	3M™ <i>Listeria</i> Enrichment Broth UVM Modified BP0255500, 500g	1/case
BP0269500	3M™ Modified Letheen Broth BP0269500, 500g	12/case
BP0270500	3M™ MRS Agar BP0270500, 500g	12/case
BP0275500	3M™ MRS Broth BP0275500, 500g	1/case
BPW010	3M™ Buffered Peptone Water Broth ISO BPW010, 10kg	1/case
BPW025	3M™ Buffered Peptone Water Broth ISO BPW025, 2.5kg	1/case
BPW500	3M™ Buffered Peptone Water Broth ISO BPW500, 500g	1/case

Sample Collection, Preparation & Processing — Process

Catalog Number	Description	Quantity/Billing Unit
General Lab Equipment		
6485	3M™ Pipettor Tip 6485, 1mL, wide bore, racked	960/case
6487	3M™ Pipettor Tip 6487, 5mL, wide bore, bulk	250/case
6488	3M™ Pipettor Tip 6488, 5mL, wide bore, racked	500/case
6489	3M™ Pipettor Tip 6489, 5mL, wide bore, racked, pre-sterilized	500/case
6486	3M™ Pipettor Tip 6486, 1mL, wide bore, racked, pre-sterilized	960/case

3M™ Microbial Luminescence System

Know quicker. Release faster.

The 3M™ Microbial Luminescence System (MLS) is a rapid method for release testing of Ultra High Temperature (UHT) and Extended Shelf Life (ESL) processed products. The 3M MLS allows faster positive release of products over 48 hours sooner than traditional methods, in a semi-automated format for high throughput applications.



3M™ Microbial Luminescence System Ultra High Temperature (UHT) Beverage Screen Kit

Delivers real-time results in 27 minutes or less when testing 96 samples. Expands the application from UHT dairy products to a wide variety of UHT and ESL beverages such as fruit juices, caffeinated drinks, coconut waters, smoothies, dairy, dairy substitutes and dairy/juices mixtures.



3M™ Microbial Luminescence System Ultra High Temperature (UHT) Dairy Screen Kit

Designed to detect a broader range of microorganisms than traditional methods. Can be used for a wide variety of dairy UHT beverages.



3M™ Microbial Luminescence System

Catalog Number	Description	Quantity/Billing Unit
3M™ Microbial Luminescence Systems		
MLSII	3M™ Microbial Luminescence System MLSII	1/case
BEV600	3M™ Microbial Luminescence System Ultra High Temperature (UHT) Beverage Screen Kit BEV600	600 tests
3060	3M™ Microbial Luminescence System Mini Ultra High Temperature (UHT) Dairy Screen Kit 3060	600/case
DPQCOG3000	3M™ Microbial Luminescence System Ultra High Temperature (UHT) Dairy Screen Kit DPQCOG3000	3000/case
3M™ Microbial Luminescence System Accessories		
3004	3M™ Microbial Luminescence System ATP Reagent Control Kit 3004	50/case
3005	3M™ Microbial Luminescence System Injector Cleaning Kit 3005	1/case
3006	3M™ Microbial Luminescence System Maintenance Solution 3006	1/case
3007	3M™ Microbial Luminescence System Microwell Plate 3007	50/case
3008	3M™ Microbial Luminescence System Microwell Strip 3008	320/case
3009	3M™ Microbial Luminescence System Microwell Strip Holder 3009	1/case
ATP10	3M™ Clean-Trace™ Surface Positive Control ATP10	10/case

3M™ Time & Temperature Indicators

3M offers a variety of temperature indicators to detect rising and falling temperatures through the entire supply chain. These products can be used in an array of industries including food, pharmaceutical, medical, veterinary and chemical industry companies.



3M™ Time & Temperature Indicators

Catalog Number	Description	Quantity/Billing Unit
3M™ Freeze Watch™		
9800FW	3M™ Freeze Watch™ Indicator 9800FW, -4°C/25°F	400/case
9805FW	3M™ Freeze Watch™ Indicator 9805FW, 0°C/32°F	400/case
3M™ MonitorMark™		
9860A	3M™ MonitorMark™ Time Temperature Indicator 9860A, -15°C/5°F, Record Time 48 hours	500/case
9860B	3M™ MonitorMark™ Time Temperature Indicator 9860B, 5°C/41°F, Record Time 48 hours	5 cartons/case
9860C	3M™ MonitorMark™ Time Temperature Indicator 9860C, 10°C/50°F, Record Time 48 hours	5 packs/case
9860D	3M™ MonitorMark™ Time Temperature Indicator 9860D, 10°C/50°F, Record Time 1 week	5 packs/case
9860E	3M™ MonitorMark™ Time Temperature Indicator 9860E, 26°C/79°F, Record Time 48 hours	5 packs/case
9860H	3M™ MonitorMark™ Time Temperature Indicator 9860H, 31°C/88°F, Record Time 1 week	5 cartons/case
9861A	3M™ MonitorMark™ Time Temperature Indicator 9861A, 10°C/50°F, Record Time 2 weeks	500/case

3M™ Food Safety Hardware Service

Maintain the performance of your 3M Food Safety Hardware.



3M™ Food Safety Hardware Service

Service Name	Includes
Standard Limited Warranty	<ul style="list-style-type: none"> • Coverage for one year, starting at purchase date • Annual Preventive Maintenance and Calibration • Repairs (not due to misuse) • Loaner hardware during annual preventive maintenance and calibration and repairs
Extended Service Agreement Must be purchased within 60 days of hardware purchase	<ul style="list-style-type: none"> • Two additional years of free annual preventive maintenance and repairs, as described in the standard warranty (3 years total)
Annual Preventive Maintenance and Calibration*	<ul style="list-style-type: none"> • One-time maintenance and calibration of your hardware • Inspection of the instrument • Cleaning • Preventive Maintenance Service • Run software functionality tests for instrument and performance • Perform annual calibration against National Institute of Standards and Technology (NIST) standard • Documentation • Loaner hardware*
Repair*	<ul style="list-style-type: none"> • One-time repair, including parts replacement and labor • Cleaning and calibration • Loaner hardware*

Terms and prices subject to change.

*Services included in warranty and extended service agreement; charges for units out of warranty.

Catalog Number	Description	Quantity/Billing Unit
Extended Service Agreements		
6558-ESA	PPRA 2YR Extended Service Agreement	1/case
MDSSERV	MDS and MLS 2YR Extended Service Agreement	1/case
UNGESA	LM1 2YR Extended Service Agreement	1/case



3M Food Safety

3M Center, Building 275-5W-05
St. Paul, MN 55144-1000 USA

1-800-328-6553
[3M.com/FoodSafety](https://www.3m.com/FoodSafety)

3M, Clean-Trace, Freeze Watch, MonitorMark and Petrifilm are trademarks of 3M.
All other trademarks are property of their respective owners. Please recycle. Printed in USA.
© 3M 2021. All rights reserved. 70-2011-5205-8 Rev-0621