



# *Lab Vortex Mixers Accessories*

*Catalog*



***Explore the various accessories  
that you can use with your vortex mixer***  
 *HOLDERS, MICROPLATES, FOAM INSERTS, AND MORE!*

***Ingeniously Practical***

***WWW.OHAUS.COM***

# VORTEX MIXER ACCESSORIES

## Micro-Tube Holder

Mixes (48) 0.25 to 2 mL micro-tubes.  
Requires Insert Retainer.

Description	Item No.
Micro-Tube Holder (2 pack)	30400232
Insert Retainer	30400227

## 9 to 13 mm Tube Holder

Ideal for mixing 5 mL culture tubes and micro-vials.  
Requires Insert Retainer.

Description	Item No.
9 to 13 mm Tube Holder (2 pack)	30400229
Insert Retainer	30400227

## 14 to 19 mm Tube Holder

Ideal for mixing up to (8) 15 mL centrifuge tubes.  
Requires Insert Retainer.

Description	Item No.
14 to 19 mm Tube Holder (2 pack)	30400230
Insert Retainer	30400227

## 20 to 25 mm Tube Holder

Ideal for mixing up to (8) 50 mL centrifuge tubes.  
Requires Insert Retainer

Description	Item No.
20 to 25 mm Tube Holder (2 pack)	30400231
Insert Retainer	30400227

## Vessel Harness

Mixes Erlenmeyer flasks and media bottles.  
Requires Insert Retainer.

Description	Item No.
Vessel Harness (2 pack)	30400228
Insert Retainer	30400227

## Cup Head

Designed for mixing 1 tube at a time.

Description	Item No.
Cup Head	30400235

## 3" Rubber Head Cover and 3" Head

Designed for mixing irregular shaped objects.

Description	Item No.
3" (7.6 cm) Rubber Head Cover	30400196
3" (7.6 cm) Head	30400236
3" (7.6 cm) Head with Rubber Cover	30400237

## Microplate Holder

Ideal for mixing 96-well plates or deep well blocks.  
Requires Insert Retainer.

Description	Item No.
Microplate Holder (2 pack)	30400233
Insert Retainer	30400227

## Flat Foam Insert

Ideal for custom applications. Can be cut or drilled to fit your specifications. Requires Insert Retainer.

Description	Item No.
Flat Foam Insert (2 pack)	30400234
Insert Retainer	30400227

## Microwell Kit with Retainer

Includes:  
• 2 Microplate Holders  
• 1 Insert Retainer

Description	Item No.
Microwell Kit with Retainer	30400199

## Microtube Kit with /Retainer

Includes:  
• 2 Micro-Tube Holders  
• 1 Insert Retainer

Description	Item No.
Microtube Kit with /Retainer	30400200

## Foam Insert Variety Pack

Includes:  
• (1) 9 to 13 mm Tube Holder  
• (1) 14 to 19 mm Tube Holder  
• (1) 20 to 25 mm Tube Holder  
• 1 Flat Foam Insert

Description	Item No.
Foam Insert Variety Pack	30400201

## Micro-Tube and Microplate Holder Kit

Includes:  
• 1 Micro-Tube Holder  
• 1 Microplate Holder  
• 1 Insert Retainer

Description	Item No.
Micro-Tube & Microplate Holder Kit	30400198

## Tube Holder Kit

Includes:  
• (1) 9 to 13 mm Tube Holder  
• (1) 14 to 19 mm Tube Holder  
• (1) 20 to 25 mm Tube Holder  
• 1 Flat Foam Insert  
• 2 Vessel Harnesses  
• 1 Insert Retainer

Description	Item No.
Tube Holder Kit	30400197

OHAUS Vortex Mixer accessories are interchangeable with other manufacturers' models.

# VORTEX MIXER ACCESSORIES

## Single Tube Holder

Single tube, hands free mixing. Easily attaches to the top of any Vortex Mixer with the use of a strong magnetic base. Accepts tubes from 2.5 to 4.5" (6.4 to 11.4 cm) in length. Minimum tube diameter of 0.75" (19 mm).

Description	Item No.
Single Tube Holder	30400202

## Single Tube Holder Accessories

(requires Single Tube Holder)

### 0.5 mL Micro-Tube Holder

Mixes (24) 0.5 mL micro-tubes.  
For use with Single Tube Holder.

Description	Item No.
Micro-Tube Holder (0.5 mL)	30400204

### Ampule Tube Holder

Mixes up to 4 storage vials and test tubes.

Description	Item No.
15 to 17 mm Ampule Tube Holder	30400207
10 to 17 mm Ampule Tube Holder	30400208

### Stainless Steel Microtube Holder

Mixes up to (12) 1.5 to 5 mL tubes.  
Stainless steel construction.

Description	Item No.
Microtube Holder	30400206

## Adapter for Vortex-Genie® Mixer\*

Adapter plate easily adheres to the Vortex-Genie® Mixer housing so Single Tube Holder (sold separately) can be attached.

Description	Item No.
Adapter plate	30400203

\*The Vortex-Genie® Mixer is a registered trademark of Scientific Industries, Inc.

### 1.5 mL to 2.0 mL Micro-Tube Holder

Mixes (18) 1.5 to 2 mL micro-tubes.  
For use with Single Tube Holder.

Description	Item No.
Microtube Holder (1.5 mL to 2.0 mL)	30400205

# HEAVY DUTY AND MICROPLATE VORTEX MIXER ACCESSORIES

## Foam Insert for 0.5 mL Microtubes

Foam insert holds (52) 0.5 mL microtubes.  
Requires Universal Holder.

Description	Item No.
0.5 mL microtubes	30400216
Universal Holder	30400226

## Foam Insert for 1.5 to 2.0 mL

### Microtubes

Foam insert holds (38) 1.5 to 2.0 mL microtubes.  
Requires Universal Holder.

Description	Item No.
1.5 to 2.0 mL microtubes	30400217
Universal Holder	30400226

## Foam Insert for 12-13 mm Test Tubes

Foam insert holds (34) 12-13 mm diameter test tubes.  
Requires Universal Holder.

Description	Item No.
12-13 mm test tubes	30400220
Universal Holder	30400226

## Foam Insert for 15-18 mm Test Tubes

Foam insert holds (20) 15-18 mm diameter test tubes.  
Ideal for 15 mL centrifuge tubes.  
Requires Universal Holder.

Description	Item No.
15-18 mm test tubes	30400221
Universal Holder	30400226

## Foam Insert for 19-21 mm Test Tubes

Foam insert holds (18) 19-21 mm diameter test tubes.  
Requires Universal Holder.

Description	Item No.
19-21 mm test tubes	30400222
Universal Holder	30400226

## Foam Insert for 22-25 mm Test Tubes

Foam insert holds (13) 22-25 mm diameter test tubes.  
Requires Universal Holder.

Description	Item No.
22-25 mm test tubes	30400223
Universal Holder	30400226

## Foam Insert for 26-29 mm Test Tubes

Foam insert holds (4) 26-29 mm diameter test tubes. Ideal for 50 mL centrifuge tubes. Requires Universal Holder.

Description	Item No.
26-29 mm test tubes	30400224
Universal Holder	30400226

## Cup Head\*

Designed for mixing 1 tube at a time.

Description	Item No.
Cup Head	30400210

## Small Vessel Holder

Rubber holder secures 125 and 250 mL Erlenmeyer flasks.  
Vessel holder also includes a grip mat.  
Requires Universal Holder.

Description	Item No.
Small Vessel Holder	30400218
Universal Holder	30400226

## Large Vessel Holder

Rubber holder secures 500 and 1000 mL Erlenmeyer flasks.  
Vessel holder also includes a grip mat.  
Requires Universal Holder.

Description	Item No.
Large Vessel Holder	30400211
Universal Holder	30400226

## Microplate Holder (Single)\*

Designed to hold one standard microplate.

Description	Item No.
Single Holder	30400215

## Microplate Holder (Double)\*

Designed to hold two standard microplates.

Description	Item No.
Double Holder	30400213

## Microplate Holder (Quad)

Designed to hold four standard microplates.

Description	Item No.
Quad Holder	30400214

## Stackable Microplate Holder Four

Designed to maximize the capacity of the Heavy-Duty Vortex Mixer to eight microplates by stacking the tray on top of the four plate holder.

Description	Item No.
Stackable Holder	30400212

## Flat Foam Insert

Ideal for custom applications. Can be cut or drilled to fit your specifications. Requires Universal Holder.

Description	Item No.
Flat Foam Insert	30400209
Universal Holder	30400226

## Universal Holder & Cover


Replacement for items supplied with Heavy-Duty Vortex Mixer. Cover allows for mixing irregularly shaped objects.

Description	Item No.
Universal Holder	30400226
Universal Holder Cover	30400225

\*Unless noted with an asterisk (\*), for use on Heavy-Duty Vortex only.

# MULTI-TUBE VORTEXER ACCESSORIES

## Foam Test Tube Racks



Note: Foam colors are subject to change

Description	Tube Capacity	Dimensions (L x W x H)	Item No.
10 mm Test Tube Foam Rack	50	5.5 x 9.5 x 2" (14 x 24.1 x 5.1 cm)	30400238
12 mm Test Tube Foam Rack	50	5.5 x 9.5 x 2" (14 x 24.1 x 5.1 cm)	30400239
13 mm Test Tube Foam Rack	50	5.5 x 9.5 x 2" (14 x 24.1 x 5.1 cm)	30400240
16 mm Test Tube Foam Rack (for 15 mL centrifuge tubes)	50	5.5 x 9.5 x 2" (14 x 24.1 x 5.1 cm)	30400241
25 mm Test Tube Foam Rack	28	5.5 x 9.5 x 2" (14 x 24.1 x 5.1 cm)	30400243
29 mm Test Tube Foam Rack (for 50 mL centrifuge tubes)	15	5.5 x 9.5 x 2" (14 x 24.1 x 5.1 cm)	30400242
Replacement Tray Pad Set (upper & lower)	N/A	7 x 12 x 1" (17.8 x 30.5 x 2.5 cm)	30400245

## Post Extension Kit

Adds 6" (15.2 cm) to the post of the Multi-Tube Vortex Mixer to accommodate tubes up to 10" (25.4 cm) tall.


Description	Item No.
Post Extension Kit	30400244





**Visit the Lab Vortex Accessories Product Page:**

Click the image or use the QR code



**Contact Us:**

**OHAUS CORPORATION**

Tel: 800.672.7722  
973.377.9000  
Fax: 973.944.7177

**OHAUS Europe GmbH**

e-mail: [ssc@ohaus.com](mailto:ssc@ohaus.com)  
Tel: 0041 22 567 53 19  
e-mail: [tsc@ohaus.com](mailto:tsc@ohaus.com)  
Tel: 0041 22 567 53 20

**OHAUS Asia/Pacific Headquarters**

e-mail: [APmarketing@ohaus.com](mailto:APmarketing@ohaus.com)  
Tel: +86 21 64855408

**奥豪斯仪器 OHAUS Instruments**

客户咨询邮箱: [inquirychina@ohaus.com](mailto:inquirychina@ohaus.com)  
售后服务邮箱: [ohauservice@ohaus.com](mailto:ohauservice@ohaus.com)  
热线电话: 4008-217-188

[www.ohaus.com](http://www.ohaus.com)

80776906\_B 20240618 © Copyright OHAUS Corporation

*The management system governing the manufacture of this product is ISO 9001:2015 certified.*

