

## Description

The SMRTbell enzyme cleanup kit 1.0 contains 4 enzymes and buffer for enzymatic removal of incomplete or damaged SMRTbell molecules (18 samples).

## Kit configuration & layout

### Configuration

#	Component	Part number	Qty	Color	Volume
1	10X clean up buffer	101-741-600	1	white	1000 µL
2	Enzyme A	101-741-100	1	yellow	80 µL
3	Enzyme B	101-741-700	1	blue	25 µL
4	Enzyme C	101-741-400	1	red	25 µL
5	Enzyme D	101-741-500	1	orange	40 µL

### Layout



## Quantities

18 reactions

## Storage and handling requirements

Product should be stored frozen between -15 to -25 °C.

Refer to protocols for additional recommendations for reagent handling.

## Shelf life

Refer to expiration date on kit label.

## Intended use

Research use only. Not for use in diagnostic procedures.

## Support

For additional product and technical information, such as safety data sheets (SDS) or procedures, visit our [website](#).

## Legal

Information in this document is subject to change without notice. PacBio assumes no responsibility for any errors or omissions in this document. Certain notices, terms, conditions and/or use restrictions may pertain to your use of PacBio products and/or third-party products. Refer to the applicable PacBio terms and conditions of sale and to the applicable license terms at [pacb.com/license](http://pacb.com/license). Pacific Biosciences, the PacBio logo, PacBio, Circulomics, Omniome, SMRT, SMRTbell, Iso-Seq, Sequel, Nanobind, and SBB are trademarks of PacBio. All other trademarks are the sole property of their respective owners.